

Making a
Difference
in the
Journey
of Life

Ostomy Care Product Catalog



Ostomy Care
Healthy skin. Positive outcomes.

Healthy Skin. Positive Outcomes.

This is our quest. Keeping peristomal skin healthy. The way it was meant to be.

Hollister Ostomy Care is dedicated to delivering products that provide an optimum fit to prevent leakage, along with the best skin interface formulations to protect and maintain healthy skin.

We believe every unique person with a stoma deserves a customized solution, and the positive outcomes that can result. Our products and services are designed to do just that.

Thank you for partnering with us on our journey.

Table of Contents

Why Secure Start Services for Your Patients.....	4-5
Choosing the Right Skin Barrier.....	8-11
Ceramide-Infused Products.....	12-23
New Image Two-Piece Pouching System	15-16
CeraPlus Flat Skin Barriers with Integrated Floating Flange.....	15
CeraPlus Convex Skin Barriers with Integrated Floating Flange...15-16	
Premier One-Piece Pouching System	17-20
Drainable Pouches with CeraPlus Flat Barriers	17
Drainable Pouches with CeraPlus Soft Convex Barriers.....	18
Drainable Pouches with CeraPlus Convex Barriers	19
Urostomy Pouches with CeraPlus Soft Convex Barriers.....	20
Non Sterile Kits - Drainable	21
Non Sterile Kits - Urostomy	21
Ostomy Accessories	23
Adapt CeraRing Flat Barrier Rings	23
Adapt CeraRing Convex Barrier Rings	23
Adapt CeraRing Oval Convex Barrier Rings	23
New Image Two-Piece Pouching System	26-46
Flat Skin Barriers with Integrated Floating Flange.....	29-31
Convex Skin Barriers with Integrated Floating Flange	32-35
Soft Convex Skin Barriers with Integrated Floating Flange.....	36
Drainable Pouches.....	38-39
High Output Drainable Pouches.....	40
Closed Pouches.....	41
Urostomy Pouches.....	42
Non Sterile Kits — Drainable	43-44
Non Sterile Kits — Urostomy	44
Sterile Kits — Drainable	45
Sterile Kits — Urostomy	46
Premier One-Piece Pouching System	47-71
Drainable Pouches with Flat Barriers	48-53
Drainable Pouches with Soft Convex Barriers.....	55-56
Drainable Pouches with Convex Barriers.....	57-62
High Output Drainable Pouches.....	63
Closed Pouches with Flat Barriers	64-65
Stoma Cap	65
Urostomy Pouches with Flat Barriers	66-67
Urostomy Pouches with Convex Barriers.....	68-69
Urostomy Pouches with CeraPlus Soft Convex Barriers.....	70
Non Sterile Kits — Drainable	71
Non Sterile Kits — Urostomy	71
Pouchkins Pediatric Pouching System	73-76
Premie Pouch.....	73
Newborn Pouch.....	73
One-Piece Pouches	73-74
Two-Piece Skin Barrier.....	75
Two-Piece Pouches	75-76
Pediatric Ostomy Belt.....	76

Ostomy Accessories	78-89
Adapt Skin Barrier Rings	78
Adapt CeraRing Barrier Rings	80
Adapt Skin Barrier Strips	80
Adapt Barrier Extenders	80
Adapt Skin Barrier Paste	80
Adapt Stoma Powder	80
Adapt Odor Eliminator and Lubricant	80
Adapt Ostomy Belts	80
Adapt Skin Prepping Product.....	81
Adapt Adhesives and Removers	81
Flextend (extended wear) Skin Barrier	81
Odor Eliminators	82
Cleaner/Decrystallizer	82
Drainable Pouch Clamp	82
Urostomy Drain Tube Adapter	82
Hollister Universal Catheter Access Port.....	82
Loop Ostomy Bridge	83
Bedside Drainage Collection System with Anti-Reflux Valve	83
Skin Barrier Pastes and Powders.....	83
Skin Barriers.....	83
Replacement Filters	83
Vented Urinary Leg Bag Combination Pack	84
Urinary Leg Bag Combination Pack	84
Urinary Leg Bag.....	84
Fabric Leg Bag Straps	85
Extension Tubing.....	85
Irrigation Products.....	86-87
Karaya 5 One-Piece Pouching System	88-89
Drainable Pouches.....	88-89
Closed Pouches.....	89
Urostomy Pouches.....	89
Education	91-97
Resources.....	92-93
Glossary	94-96
Enterocutaneous Fistula Management	97
Stock Number Index	98-99
Product Name Index	100-101
Reimbursement for Selected Ostomy Products	102-103
Medicare Utilization Chart	104-107
HCPCS Codes Index	108

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Hollister Secure Start Services

Personalized Patient Support



Support Throughout the Continuum of Care

We are proud to offer dedicated support for each and every patient we serve along the continuum of care. Hollister Secure Start services provide a lifetime of personalized support for our members.

Our team is focused on helping patients with:

- Finding the right product fit
- Navigating insurance coverage
- Identifying supplier options
- Providing condition specific education
- And more!

To learn more, call us at **1.888.808.7456**.



The Hollister logo and Secure Start are trademarks of Hollister Incorporated.
© 2020 Hollister Incorporated.

Why Secure Start Services for Your Patients?

Secure Start services helps patients—regardless of product brand they use—to get back to living life with confidence and security. We understand this is a process that takes time. That's why our no-cost, personalized service is available for the lifetime of their ostomy.



THE RIGHT PRODUCTS—based on your product recommendations, we help patients find the ostomy products that best fit their needs throughout the lifetime of their ostomy.



OSTOMY SMART EDUCATION—easy access to important care information and lifestyle resources that support the goal of maintaining healthy peristomal skin from the start.



SUPPLIER MATCHING—connects patients with a supplier of their choice based on their insurance plan for the products you recommend.



WOC ACCESS—if the need arises, an ostomy nurse is available to assist with product questions.

ENROLLMENT MADE EASY

Obtaining written patient consent to enroll in Secure Start services is an important step. At Hollister, we prioritize patient privacy—giving patients comfort that their information will remain confidential and will be shared only with those they wish to hear from. Having patient consent on file helps us provide proactive ostomy education and product information throughout their journey.

Digital

Easy mobile enrollment with Enroll on the Go



Fax

Fax enrollment form to 1.800.398.2341.

If you have any questions about the enrollment process, please call your Hollister Secure Start services coordinator at 1.888.808.7456. Learn more at www.securestartservices.com.

How to Use The Hollister Product Catalog

Hollister has created this simplified catalog to help our customers find what they need in a quick and efficient manner.

If you...

know the stock number (or stock series number)

Go to the **numeric index at the back of the catalog** on pages 100-101 to find the page number where it appears.

If you...

know the product brand or series name

Go to the **Table of Contents** on pages 2-3 to find the page number where it appears or reference the Product Name Index on pages 98-99.

If you...

know the HCPCS code

Go to the **HCPCS Code Index** at the back of the catalog on page 108.

Fit and Formulation

Promoting healthy peristomal skin

It drives everything we do. In order to promote healthy skin, it is important to understand that the *fit* of the pouching system as well as the *formulation* of the skin barrier go hand in hand.

A pouching system will not maintain a good seal without the proper *fit*. Without the correct *formulation*, the skin around the stoma can be at risk for complications.

Our goal is to help keep the skin around your patient's stoma as healthy as the skin on the rest of their body.

Choosing the Right Skin Barrier

FIT is critical to get a good seal around the stoma to prevent leakage

Barrier Examples



Flat Skin Barrier



Flat Skin Barrier with Round or Oval Convex Barrier Ring

Oval barrier ring shown in picture.



Integrated Soft Convex Barrier



Integrated Convex Barrier



Integrated Convex Barrier with Round or Oval Convex Barrier Ring

Round barrier ring shown in picture.

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.

Possible Clinical Case Scenarios



Flat Skin Barrier

- Flat peristomal plane
- Adequate stoma protrusion
- Centered stoma opening



Flexible Convexity

- Stoma flush with skin
- Uneven peristomal plane
- Off-centered stoma opening



Soft Convexity

- Gentle pressure required to peristomal plane
- Flush stoma
- Peristomal injury



Firm Convexity

- Creases or folds in peristomal plane
- Off-centered stoma opening
- Stoma retracted below skin surface



Firm Convexity with More Depth

- Recessed or flushed stoma
- Uneven abdominal contours
- Off-centered stoma opening

Choosing the Right Skin Barrier Formulation

FORMULATION is the customized ingredients responsible for adhesion, erosion resistance and absorption that must be carefully formulated to create an environment where healthy skin can thrive.

Types of Skin Barrier Formulations and Their Performance

The Hollister Ostomy Care portfolio has a multitude of formulation options that are designed to help ensure an optimal skin interface is possible for every patient.

Extended Wear

CeraPlus Skin Barrier

- Contains all skin barrier characteristics (adhesion, absorption and erosion resistance) and is infused with ceramide – a naturally occurring lipid in the skin which protects against dryness
- Its formulation is designed to help protect the skin's natural moisture barrier and help maintain good peristomal skin health

Flexend Skin Barrier

- Durable with high resistance to erosion
- Swells when it comes in contact with liquid to help form a seal between the skin and the stoma
- High initial tack and suitable for those who perspire heavily
- Strong adhesion and flexibility

FormaFlex Skin Barrier

- Helps to provide a secure peristomal skin seal as it will maintain its shape after being stretched or shaped
- Provides an easy-to-customize opening due to barrier material
- Easy to shape to individual stomas without the use of scissors

Standard Wear

FlexWear Skin Barrier

- Good balance of durability and gentleness providing moderate adhesion, and can handle frequent or infrequent pouch changes
- Resistant to fluid erosion and tends to hold its shape in the presence of fluids
- Suitable for all skin types, depending on change frequency
- Designed to be worn more than one day

SoftFlex Skin Barrier

- Appropriate choice for:
 - Those who have thin, fragile, or sensitive skin (such as those who have undergone chemo/radiation therapy)
 - People who desire frequent removal
- Designed to limit skin damage as its gentle to the skin, which allows for frequent pouch changes
- Thin and flexible without compromising security

Karaya 5 Skin Barrier

- Well-suited for sensitive skin or in some cases where allergic reactions to synthetic barriers has occurred
- Mildly adhesive – needs a belt to secure pouching system

Combined with the flexible and customizable fit choices we offer in our pouching systems, barrier styles and accessories, you can count on the right combination of fit and formulation for your patients.

Pouch and barrier fit options come in both 1 piece and 2 piece options.

- Flat
- Convex
- Soft Convex

The CeraPlus skin barrier

helps maintain healthy skin

from Day One



Ceramide Products

Transforming the skin health conversation


At Hollister Ostomy Care, we believe peristomal skin deserves better. In fact, we think we should expect the skin around the stoma to be as healthy as the skin on the other side of the abdomen.

That's why we've made the pursuit of healthy skin our priority.

The CeraPlus skin barrier is formulated to support healthy skin from Day One. The ceramide-infused barrier is designed to help protect the skin's natural moisture barrier. Healthy skin can remain healthy through pouching changes, and support an appropriate fit around the stoma.

Maintaining the Right Environment

Makes all the Difference



**Available on Premier One-Piece and
New Image Two-Piece Pouching Systems**

The CeraPlus Skin Barrier with Remois Technology*

At Hollister, we realize the importance of healthy peristomal skin. Our newest skin barrier is infused with ceramide, the skin's naturally occurring protection against dryness. The CeraPlus skin barrier is designed to help protect the skin's own moisture barrier, and is available on Premier one-piece and New Image two-piece pouching systems. This means you can enjoy the same pouch and skin barrier features from Hollister that you've come to rely on—and more!

The CeraPlus skin barrier from Hollister Ostomy Care—because peristomal skin deserves better.



*Remois is a technology of Alcare Co., Ltd.



Ostomy Care
Healthy skin. Positive outcomes.

Ceramide-Infused Products

The CeraPlus Skin Barrier delivers both Fit and Formulation

Ceramide-Infused Barrier

- Protects the skin's natural moisture barrier
- Decreases TEWL (transepidermal water loss) from damaged or eroded skin
- Maintains peristomal skin health

Color Matching System

Two-Piece Pouching Systems

To help you select compatible *New Image* pouching system components, every box label is color-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its color code.
2. Select a suitable pouch with the desired features in the matching color.

Flange Sizes:



1-3/4" (44 mm)



2-1/4" (57 mm)



2-3/4" (70 mm)



4" (102 mm)

NOTE: All stock numbers shaded in gray throughout this catalog denote ceramide-infused products.

Look for the convex and soft convex symbols:



CeraPlus Flat Skin Barriers with Integrated Floating Flange



1120x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5
Ceramide

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
11202	1-3/4" (44 mm)	1-1/4" (32 mm)
11203	2-1/4" (57 mm)	1-3/4" (44 mm)
11204	2-3/4" (70 mm)	2-1/4" (57 mm)



1510x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5
Ceramide

HCPCS for 15102, 15103: A4409

HCPCS for 15104: A4410

Stock No	Flange Size	For Stomas Up To
15102	1-3/4" (44 mm)	1-1/4" (32 mm)
15103	2-1/4" (57 mm)	1-3/4" (44 mm)
15104	2-3/4" (70 mm)	2-1/4" (57 mm)

CeraPlus Convex Skin Barriers with Integrated Floating Flange



1140x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border

CONVEX
Ceramide

HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11402	1-3/4" (44 mm)	1" (25 mm)
11403	2-1/4" (57 mm)	1-1/2" (38 mm)
11404	2-3/4" (70 mm)	2" (51 mm)



1530x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Without Tape Border

CONVEX
Ceramide

HCPCS for 15302, 15303: A4407

HCPCS for 15304: A4408

Stock No	Flange Size	For Stomas Up To
15302	1-3/4" (44 mm)	1" (25 mm)
15303	2-1/4" (57 mm)	1-1/2" (38 mm)
15304	2-3/4" (70 mm)	2" (51 mm)



1150x (PRE-SIZED)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ With Tape Border

CONVEX
Ceramide

HCPCS: A4407

Stock No	Flange Size	Barrier Opening
11503	1-3/4" (44 mm)	7/8" (22 mm)
11504	1-3/4" (44 mm)	1" (25 mm)
11505	2-1/4" (57 mm)	1-1/8" (29 mm)
11506	2-1/4" (57 mm)	1-1/4" (32 mm)

Ceramide-Infused Products

CeraPlus Convex Skin Barriers with Integrated Floating Flange (continued)



New!

1170x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Soft Convex
- ▶ With Tape Border

Box of 5



Ceramide

HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11702	 1-3/4" (44 mm)	1" (25 mm)
11703	 2-1/4" (57 mm)	1-1/2" (38 mm)
11704	 2-3/4" (70 mm)	2" (51 mm)



New!

1190x (PRE-SIZED)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Soft Convex
- ▶ With Tape Border

Box of 5



Ceramide

HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11903	 1-3/4" (44 mm)	7/8" (22 mm)
11904	 1-3/4" (44 mm)	1" (25 mm)
11905	 2-1/4" (57 mm)	1-1/8" (29 mm)
11906	 2-1/4" (57 mm)	1-1/4" (32 mm)

The CeraPlus Skin Barrier delivers both Fit and Formulation

Ceramide-Infused Barrier

- **Protects** the skin's natural moisture barrier
- **Decreases** TEWL (Transepidermal water loss) from damaged or eroded skin
- **Maintains** peristomal skin health



Premier One-Piece Pouching System

Drainable Pouches with CeraPlus Flat Barriers



8931 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 10
Ceramide

HCPCS: A4388

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

8931	2-1/2" (64 mm)	12" (30 cm)
------	----------------	-------------



8901 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

Box of 10
Ceramide

HCPCS: A5056

Beige w/Viewing Option Stock No

For Stomas Up To

Pouch Length

8901	2-1/2" (64 mm)	12" (30 cm)
------	----------------	-------------



89xx (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

Box of 10
Ceramide

HCPCS: A5056

Beige w/Viewing Option Stock No

Barrier Opening

Pouch Length

8925	1" (25 mm)	12" (30 cm)
8930	1-3/16" (30 mm)	12" (30 cm)
8935	1-3/8" (35 mm)	12" (30 cm)
8940	1-9/16" (40 mm)	12" (30 cm)

Ceramide-Infused Products

Drainable Pouches with **CeraPlus** Soft Convex Barriers



895xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8954	1" (25 mm)	12" (30 cm)
8958	1-1/2" (38 mm)	12" (30 cm)
89511	2-1/8" (55 mm)	12" (30 cm)



896x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8960	3/4" (19 mm)	12" (30 cm)
8961	7/8" (22 mm)	12" (30 cm)
8962	1" (25 mm)	12" (30 cm)
8963	1-1/8" (29 mm)	12" (30 cm)

Drainable Pouches with CeraPlus Convex Barriers



891xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8914	1" (25 mm)	12" (30 cm)
8918	1-1/2" (38 mm)	12" (30 cm)
89111	2" (51 mm)	12" (30 cm)



898xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5



HCPCS Code: A4390

Ultra-Clear

Stock No	For Stomas Up To	Pouch Length
8984	1" (25 mm)	12" (30 cm)
8988	1-1/2" (38 mm)	12" (30 cm)
89811	2" (51 mm)	12" (30 cm)



899x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8990	3/4" (19 mm)	12" (30 cm)
8991	7/8" (22 mm)	12" (30 cm)
8992	1" (25 mm)	12" (30 cm)
8993	1-1/8" (29 mm)	12" (30 cm)

Ceramide-Infused Products

Urostomy Pouches with **CeraPlus** Soft Convex Barriers

Upgraded Design



8413xx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
includes 1 adapter

**SOFT
CONVEX**
Ceramide

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Ultra Clear Stock No	For Stomas Up To	Pouch Length
84134	1" (25 mm)	9" (23 cm)
84138	1-1/2" (38 mm)	9" (23 cm)
841311	2-1/8" (55 mm)	9" (23 cm)

Upgraded Design



841x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
includes 1 adapter

**SOFT
CONVEX**
Ceramide

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4430

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8412	3/4" (19 mm)	9" (23 cm)
8413	7/8" (22 mm)	9" (23 cm)
8414	1" (25 mm)	9" (23 cm)
8415	1-1/8" (29 mm)	9" (23 cm)

Non Sterile Kits — Drainable



1960x




Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear Pouch Film
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

Box of 5

Ceramide

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19602	 1-3/4" (44 mm)	1-1/4" (32 mm)
19603	 2-1/4" (57 mm)	1-3/4" (44 mm)
19604	 2-3/4" (70 mm)	2-1/4" (57 mm)



1980x




Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear Pouch Film
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only
- ▶ Clamp Included

Box of 5

Ceramide

HCPCS: A5063 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19802	 1-3/4" (44 mm)	1-1/4" (32 mm)
19803	 2-1/4" (57 mm)	1-3/4" (44 mm)
19804	 2-3/4" (70 mm)	2-1/4" (57 mm)

Non Sterile Kits — Urostomy



1990x




Urostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1892x with *ComfortWear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

Box of 5

Ceramide

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19902	 1-3/4" (44 mm)	1-1/4" (32 mm)
19903	 2-1/4" (57 mm)	1-3/4" (44 mm)
19904	 2-3/4" (70 mm)	2-1/4" (57 mm)

Infused with the Science

of Skin Health



Adapt CeraRing Barrier Ring

Sealed with the confidence of ceramide

Adapt barrier ring technology was introduced more than 10 years ago. Its ability to be stretched and shaped to help minimize the risk of leakage and protect the peristomal skin makes it a favorite among WOC nurses.



Now, innovation has reached the next level, with the introduction of Adapt CeraRing.

Infused with ceramide—a natural component of the skin—Adapt CeraRing is designed to support the healthy peristomal skin your patients deserve. Ceramides help protect the skin's natural moisture barrier and help prevent water loss that may lead to damage and dryness.

Experience the confidence of ceramide—try Adapt CeraRing today.



Ostomy Care
Healthy skin. Positive outcomes.

Ceramide-Infused Products

Ostomy Accessories

Adapt CeraRing Flat Barrier Rings



88xx

- ▶ Barrier Rings (Flat)
- ▶ *Flexlend M* (extended wear) with Ceramide

Box of 10

Ceramide

HCPCS: A4385

Stock No	Outer Diameter	Thickness
8805	2" (48 mm)	4.5 mm
8815	2" (48 mm)	2.3 mm SLIM

Adapt CeraRing Convex Barrier Rings



895xx

- ▶ Round Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10

CONVEX
Ceramide

HCPCS: A4411

Stock No	Inner Diameter
89520	13/16" (20 mm) — can be stretched to 1" (25 mm)
89530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
89540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)

Adapt CeraRing Oval Convex Barrier Rings



896xx

- ▶ Oval Convex Barrier Rings
- ▶ *Flexlend* (extended wear) with Ceramide

Box of 10

CONVEX
Ceramide

HCPCS: A4411

Stock No	Inner Diameter
89601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
89602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
89603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.



Follow the evidence to help ensure healthy skin

After surgery: Peristomal skin starts healthy



Top layers of the epidermis are removed
with *every barrier change*.

It's unavoidable, and the damage isn't
always visible.




Epidermal stripping is linked to increased
levels of transepidermal water loss
(TEWL)^{1,2}

References:

1. McNichol L, Lund C, Rosen T, Gray M. Medical adhesives and patient safety: state of the science: consensus statements for the assessment, prevention, and treatment of adhesive-related skin injuries. *J Wound Ostomy Continence Nurs.* 2013;40(4):365-380.
2. Grove G, Houser T, Sibbald G, Salvadalena G. Measuring epidermal effects of ostomy skin barriers. *Skin Res Technol.* 2019;25:179-186. <https://doi.org/10.1111/srt.12630>
3. Nybaek H, Lophagen S, Karlsmark T, Bang Knudsen D, Jemec GBE. Stratum corneum integrity as a predictor for peristomal skin problems in ostomates. *Br J Dermatol.* 2010;162:357-361.

Healthy Skin: Keeps moisture in and irritants out




Low TEWL
value = higher
skin integrity³

Ceramides are abundant in healthy skin.
They form a waterproof barrier to help keep
moisture in and bacteria out.

— = Ceramides

Damaged Skin: Allows moisture out and irritants in



Higher TEWL
value = lower
skin integrity³

Damaged skin contains fewer ceramides.

💧 = Water droplet illustrates moisture loss from skin.

Clinical Trial Results

The Results are Clear: CeraPlus Skin Barrier Helps Ensure Healthy Skin From Day One

At Hollister Ostomy Care, we believe peristomal skin deserves better. Clinical evidence continues to empower everything we do.

The ADVOCATE clinical trial looked at difference in ostomy skin barriers and their impact on costs of stoma care and occurrence of peristomal skin complications (PSCs).

- 14% relative reduction in stoma-related cost of care with the CeraPlus skin barrier

See the Evidence That Can Empower You.



CeraPlus skin barriers help support healthy skin and help minimize ostomy burden right from the start.¹



CeraPlus skin barriers significantly reduce the stoma-related cost of care.¹



Significantly more participants were "very satisfied" with overall barrier performance and itching prevention.¹



15% fewer PSCs with the CeraPlus Skin Barrier.



	CeraPlus Skin Barrier	SenSura Mio
Observed to be better at helping maintain intact skin. ^{2,3}		
Observed to be better at managing TEWL. ^{2,3}		
Observed to be better at protecting skin against damage from repeated barrier removal. ^{2,3}		

Start using the CeraPlus skin barrier from day one to help ensure the healthy skin your patients deserve.

Contact your Hollister sales representative today for details.



References:

1. Colwell J, Pittman J, Salvadelen G. A randomized controlled trial determining variances in ostomy skin conditions and the economic impact (ADVOCATE). *Journal Wound Ostomy Continence Nurs.* 2018;45(1):37-42
2. Nichols T, Houser T, Grove G. Comparing the skin stripping effects of three ostomy skin barriers infused with ceramid, honey or aloe. *Journal of Stomal Therapy Australia* 2019; 39(2); 14-18
3. Grove G, Houser T, Sibbald G, Salvadelen G. Measuring epidermal effects of ostomy skin barriers. *Skin Res Technol.* 2019;25:179-186. <https://doi.org/10/1111/srt/12630>

Keeping peristomal skin

healthy is the inspiration

for all we do

New Image
Two-Piece Pouching System



New Image
Two-Piece Pouching System

Why choose a two-piece pouching system?

With a two-piece pouching system, the skin barrier and pouch are separate and connect together with a plastic ring, or flange. Here are a few reasons why you might want to use a two-piece pouching system:

- You may want to change your pouch more often than your skin barrier
- You may want to alternate between different pouch styles. For example, you might want to wear a smaller pouch during exercise or in the shower, and a larger pouch at night
- You may find it easier to center your skin barrier over your stoma using a two-piece system



Integrated Floating Flange

- Audible clicks signal a secure attachment
- Discreet low profile attachment of pouch to barrier



CeraPlus Skin Barrier

The CeraPlus skin barrier with Remois technology* is infused with ceramide, the skin's naturally occurring protection against dryness. It is designed to maintain adhesive properties, and features a ceramide-infused formulation to help:

- Protect the skin's natural moisture barrier
- Maintain good peristomal skin health
- Decrease transepidermal water loss (TEWL) from damaged or eroded skin

**Remois technology is a technology of Alcare Co., Ltd.*

New Image Pouches



- Waterproof, integrated AF300 filter is designed to allow gas and air to escape, but keeps odor and output in the pouch. Liquid cannot penetrate the filter, therefore not requiring a protective sticker
- Choice of secure pouch closure options – clamp or integrated closure
- ComfortWear panels provide enhanced discretion and comfort
- Numerous pouch size options to meet individual preferences

New Image Two-Piece Pouching System

Color Matching System

Two-Piece Pouching Systems

To help you select compatible New Image pouching system components, every box label is color-coded.

1. Select the skin barrier flange size with the appropriate barrier stoma opening and note its color code.
2. Select a suitable pouch with the desired features in the matching color.

Flange Sizes:



1-3/4" (44 mm)



2-1/4" (57 mm)



2-3/4" (70 mm)



4" (102 mm)

NOTE: All stock numbers shaded in gray throughout this catalog denote ceramide-infused products.

Look for the convex and soft convex symbols:



For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Flat Skin Barriers with Integrated Floating Flange



1120x (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5

Ceramide

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
11202	1-3/4" (44 mm)	1-1/4" (32 mm)
11203	2-1/4" (57 mm)	1-3/4" (44 mm)
11204	2-3/4" (70 mm)	2-1/4" (57 mm)



1510x (CUT-TO-FIT)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5

Ceramide

HCPCS for 15102, 15103: A4409

HCPCS for 15104: A4410

Stock No	Flange Size	For Stomas Up To
15102	1-3/4" (44 mm)	1-1/4" (32 mm)
15103	2-1/4" (57 mm)	1-3/4" (44 mm)
15104	2-3/4" (70 mm)	2-1/4" (57 mm)



1460x (CUT-TO-FIT)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5

HCPCS for 14602, 14603 & 14604: A4409

HCPCS for 14606: A4410

Stock No	Flange Size	For Stomas Up To
14602	1-3/4" (44 mm)	1-1/4" (32 mm)
14603	2-1/4" (57 mm)	1-3/4" (44 mm)
14604	2-3/4" (70 mm)	2-1/4" (57 mm)
14606	4" (102 mm)	3-1/2" (89 mm)



1560x (CUT-TO-FIT)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5

HCPCS for 15602 & 15603: A4409

HCPCS for 15604: A4410

Stock No	Flange Size	For Stomas Up To
15602	1-3/4" (44 mm)	1-1/4" (32 mm)
15603	2-1/4" (57 mm)	1-3/4" (44 mm)
15604	2-3/4" (70 mm)	2-1/4" (57 mm)



1420x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5

HCPCS for 14202, 14203 & 14204: A4414

HCPCS for 14206: A4415

Stock No	Flange Size	For Stomas Up To
14202	1-3/4" (44 mm)	1-1/4" (32 mm)
14203	2-1/4" (57 mm)	1-3/4" (44 mm)
14204	2-3/4" (70 mm)	2-1/4" (57 mm)
14206	4" (102 mm)	3-1/2" (89 mm)

New Image Two-Piece Pouching System

Flat Skin Barriers with Integrated Floating Flange (continued)



1520x (CUT-TO-FIT)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

HCPCS for 15202 & 15203: A4414

HCPCS for 15204: A4415

Stock No	Flange Size	For Stomas Up To
15202	1-3/4" (44 mm)	1-1/4" (32 mm)
15203	2-1/4" (57 mm)	1-3/4" (44 mm)
15204	2-3/4" (70 mm)	2-1/4" (57 mm)



1470x (PRE-SIZED)

Box of 5

- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	Barrier Opening
14701	1-3/4" (44 mm)	5/8" (16 mm)
14702	1-3/4" (44 mm)	3/4" (19 mm)
14703	1-3/4" (44 mm)	7/8" (22 mm)
14704	1-3/4" (44 mm)	1" (25 mm)
14705	1-3/4" (44 mm)	1-1/8" (29 mm)
14706	1-3/4" (44 mm)	1-1/4" (32 mm)
14707	2-1/4" (57 mm)	1-3/8" (35 mm)
14708	2-1/4" (57 mm)	1-1/2" (38 mm)
14709	2-1/4" (57 mm)	1-3/4" (44 mm)



1610x (PRE-SIZED)

Box of 5

- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ Without Tape Border

HCPCS: A4409

Stock No	Flange Size	Barrier Opening
16101	1-3/4" (44 mm)	5/8" (16 mm)
16102	1-3/4" (44 mm)	3/4" (19 mm)
16103	1-3/4" (44 mm)	7/8" (22 mm)
16104	1-3/4" (44 mm)	1" (25 mm)
16105	1-3/4" (44 mm)	1-1/8" (29 mm)
16106	1-3/4" (44 mm)	1-1/4" (32 mm)
16107	2-1/4" (57 mm)	1-3/8" (35 mm)
16108	2-1/4" (57 mm)	1-1/2" (38 mm)



DID YOU KNOW...

- Each box of **New Image** pouches includes one disposal bag per pouch for discreet disposal
- Cutting guides are provided for all cut-to-fit barriers
- It is not recommended that you cut a pre-sized skin barrier
- It is not recommended to cut past the stop line on a cut-to-fit barrier
- Available on www.hollister.com

Flat Skin Barriers with Integrated Floating Flange (continued)



1430x (PRE-SIZED)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4414

Stock No	Flange Size	Barrier Opening
14301	1-3/4" (44 mm)	5/8" (16 mm)
14302	1-3/4" (44 mm)	3/4" (19 mm)
14303	1-3/4" (44 mm)	7/8" (22 mm)
14304	1-3/4" (44 mm)	1" (25 mm)
14305	1-3/4" (44 mm)	1-1/8" (29 mm)
14306	1-3/4" (44 mm)	1-1/4" (32 mm)
14307	2-1/4" (57 mm)	1-3/8" (35 mm)
14308	2-1/4" (57 mm)	1-1/2" (38 mm)
14309	2-1/4" (57 mm)	1-3/4" (44 mm)



1640x (PRE-SIZED)

Box of 5

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

HCPCS: A4414

Stock No	Flange Size	Barrier Opening
16402	1-3/4" (44 mm)	3/4" (19 mm)
16403	1-3/4" (44 mm)	7/8" (22 mm)
16404	1-3/4" (44 mm)	1" (25 mm)
16405	1-3/4" (44 mm)	1-1/8" (29 mm)
16406	1-3/4" (44 mm)	1-1/4" (32 mm)
16407	2-1/4" (57 mm)	1-3/8" (35 mm)
16408	2-1/4" (57 mm)	1-1/2" (38 mm)



1410x (SHAPE-TO-FIT)

Box of 5

- ▶ *FormaFlex* (extended wear) Skin Barrier
- ▶ With Tape Border

HCPCS: A4409

Stock No	Flange Size	For Stomas Up To
14102	1-3/4" (44 mm)	1/2"-1-1/4" (13-32 mm)
14103	2-1/4" (57 mm)	5/8"-1-11/16" (16-43 mm)
14104	2-3/4" (70 mm)	5/8"-2-1/4" (16-57 mm)



DID YOU KNOW...

A True Original

The Hollister integrated floating flange—as featured in the *New Image* two-piece pouching system—offers optimal flexibility and security for you and your patients. The first of its kind, this integrated floating flange helps relieve pressure on the tender postoperative or soft abdomen while offering an easy-to-use, low-profile solution patients can continue to use after they transition home. Included in all flat and convex *New Image* skin barriers.



New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange



- 1140x** (CUT-TO-FIT)
 ▶ *CeraPlus* (extended wear) Skin Barrier
 ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11402	1-3/4" (44 mm)	1" (25 mm)
11403	2-1/4" (57 mm)	1-1/2" (38 mm)
11404	2-3/4" (70 mm)	2" (51 mm)



- 1530x** (CUT-TO-FIT)
 ▶ *CeraPlus* (extended wear) Skin Barrier
 ▶ Without Tape Border

Box of 5



HCPCS for 15302, 15303: A4407

HCPCS for 15304: A4408

Stock No	Flange Size	For Stomas Up To
15302	1-3/4" (44 mm)	1" (25 mm)
15303	2-1/4" (57 mm)	1-1/2" (38 mm)
15304	2-3/4" (70 mm)	2" (51 mm)



- 1480x** (CUT-TO-FIT)
 ▶ *Flexend* (extended wear) Skin Barrier
 ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
14802	1-3/4" (44 mm)	1" (25 mm)
14803	2-1/4" (57 mm)	1-1/2" (38 mm)
14804	2-3/4" (70 mm)	2" (51 mm)



- 1580x** (CUT-TO-FIT)
 ▶ *Flexend* (extended wear) Skin Barrier
 ▶ Without Tape Border

Box of 5



HCPCS for 15802 & 15803: A4407

HCPCS for 15804: A4408

Stock No	Flange Size	For Stomas Up To
15802	1-3/4" (44 mm)	1" (25 mm)
15803	2-1/4" (57 mm)	1-1/2" (38 mm)
15804	2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

The **CeraPlus** skin barrier is infused with ceramide



- Ceramide is a natural component of the human skin
- Ceramide links the cells of the stratum corneum together to form a waterproof, protective barrier
- Ceramide helps prevent water loss that can lead to skin damage and dryness

Convex Skin Barriers with Integrated Floating Flange (continued)

Upgraded Design Recommended replacement for 1440x series below



1340x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border



Box of 5



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
13402	 1-3/4" (44 mm)	1" (25 mm)
13403	 2-1/4" (57 mm)	1-1/2" (38 mm)
13404	 2-3/4" (70 mm)	2" (51 mm)

Original Design




1440x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
14402	 1-3/4" (44 mm)	1" (25 mm)
14403	 2-1/4" (57 mm)	1-1/2" (38 mm)
14404	 2-3/4" (70 mm)	2" (51 mm)





1540x (CUT-TO-FIT)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	For Stomas Up To
15402	 1-3/4" (44 mm)	1" (25 mm)
15403	 2-1/4" (57 mm)	1-1/2" (38 mm)
15404	 2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

Select **New Image** two-piece pouching systems feature two important product upgrades:

- **Lock 'n Roll** microseal closure with ribbed strips to facilitate ease of use, and transparent flap for increased security
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - Flexible release liner covers both the skin barrier and the tape border to help provide improved ease of use
 - Tapered barrier edges to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

New Image Two-Piece Pouching System

Convex Skin Barriers with Integrated Floating Flange (continued)



1150x (PRE-SIZED)

- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
11503	1-3/4" (44 mm)	7/8" (22 mm)
11504	1-3/4" (44 mm)	1" (25 mm)
11505	2-1/4" (57 mm)	1-1/8" (29 mm)
11506	2-1/4" (57 mm)	1-1/4" (32 mm)

Upgraded Design Recommended replacement for 149xx series below



139xx (PRE-SIZED)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edge

HCPCS: A4407

Stock No	Flange Size	Barrier Opening
13901	1-3/4" (44 mm)	5/8" (16 mm)
13902	1-3/4" (44 mm)	3/4" (19 mm)
13903	1-3/4" (44 mm)	7/8" (22 mm)
13904	1-3/4" (44 mm)	1" (25 mm)
13905	2-1/4" (57 mm)	1-1/8" (29 mm)
13906	2-1/4" (57 mm)	1-1/4" (32 mm)
13907	2-1/4" (57 mm)	1-3/8" (35 mm)
13908	2-1/4" (57 mm)	1-1/2" (38 mm)
13909	2-3/4" (70 mm)	1-5/8" (41 mm)
13910	2-3/4" (70 mm)	1-3/4" (44 mm)
13911	2-3/4" (70 mm)	2" (51 mm)

Original Design

149xx (PRE-SIZED)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
14901	1-3/4" (44 mm)	5/8" (16 mm)
14902	1-3/4" (44 mm)	3/4" (19 mm)
14903	1-3/4" (44 mm)	7/8" (22 mm)
14904	1-3/4" (44 mm)	1" (25 mm)
14905	2-1/4" (57 mm)	1-1/8" (29 mm)
14906	2-1/4" (57 mm)	1-1/4" (32 mm)
14907	2-1/4" (57 mm)	1-3/8" (35 mm)
14908	2-1/4" (57 mm)	1-1/2" (38 mm)
14909	2-3/4" (70 mm)	1-5/8" (41 mm)
14910	2-3/4" (70 mm)	1-3/4" (44 mm)
14911	2-3/4" (70 mm)	2" (51 mm)



DID YOU KNOW...

Convex Skin Barrier Options

The *New Image* two-piece pouching system delivers several skin barrier options to fit a variety of patient needs. Built-in convexity options provide extra security for recessed and flush stomas, and peristomal skin with soft muscle tone. Convex skin barriers are also available with and without tape.

Convex Skin Barriers with Integrated Floating Flange (continued)



1590x (PRE-SIZED)

- ▶ *FlexTend* (extended wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS: A4407

Stock No	Flange Size	Barrier Opening
15903	1-3/4" (44 mm)	7/8" (22 mm)
15904	1-3/4" (44 mm)	1" (25 mm)
15905	2-1/4" (57 mm)	1-1/8" (29 mm)
15906	2-1/4" (57 mm)	1-1/4" (32 mm)

Upgraded Design Recommended replacement for 1450x series below



1350x (PRE-SIZED)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edge

HCPCS: A4373

Stock No	Flange Size	Barrier Opening
13503	1-3/4" (44 mm)	7/8" (22 mm)
13504	1-3/4" (44 mm)	1" (25 mm)
13505	2-1/4" (57 mm)	1-1/8" (29 mm)
13506	2-1/4" (57 mm)	1-1/4" (32 mm)

Original Design

1450x (PRE-SIZED)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	Barrier Opening
14503	1-3/4" (44 mm)	7/8" (22 mm)
14504	1-3/4" (44 mm)	1" (25 mm)
14505	2-1/4" (57 mm)	1-1/8" (29 mm)
14506	2-1/4" (57 mm)	1-1/4" (32 mm)



1550x (PRE-SIZED)

- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ Without Tape Border

Box of 5



HCPCS: A4373

Stock No	Flange Size	Barrier Opening
15503	1-3/4" (44 mm)	7/8" (22 mm)
15504	1-3/4" (44 mm)	1" (25 mm)
15505	2-1/4" (57 mm)	1-1/8" (29 mm)
15506	2-1/4" (57 mm)	1-1/4" (32 mm)

New Image Two-Piece Pouching System

Soft Convex Skin Barriers with Integrated Floating Flange



1170x (CUT-TO-FIT)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Soft Convex
- ▶ With Tape Border

Box of 5



Ceramide

New!

HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11702	1-3/4" (44 mm)	1" (25 mm)
11703	2-1/4" (57 mm)	1-1/2" (38 mm)
11704	2-3/4" (70 mm)	2" (51 mm)



1190x (PRE-SIZED)

- ▶ CeraPlus (extended wear) Skin Barrier
- ▶ Soft Convex
- ▶ With Tape Border

Box of 5



Ceramide

New!

HCPCS: A4407

Stock No	Flange Size	For Stomas Up To
11903	1-3/4" (44 mm)	7/8" (22 mm)
11904	1-3/4" (44 mm)	1" (25 mm)
11905	2-1/4" (57 mm)	1-1/8" (29mm)
11906	2-1/4" (57 mm)	1-1/4" (32 mm)

Two-Piece Soft Convex CeraPlus skin barrier



Scan Me

CeraPlus skin barrier now comes in two-piece soft convexity designed to achieve a better fit to help minimize leakage and a formulation to support healthy skin from Day One.

Watch the video to learn more about the product.



DID YOU KNOW...



Get FIT

with the digital ostomy skin barrier Fit Indicator Tool

Developed with an international panel of experienced stoma care nurses, a new tool from Hollister Education has been designed to provide a simple patient assessment methodology to determine the need for a flat or convex barrier.

Ready to get started? Visit www.hollister.com/education or download the app by searching "Hollister Education" on iTunes or Google Play.



Choose the skin barrier that

delivers Fit and Formulation to help

ensure healthy skin from Day One

The **Two-Piece Soft Convex CeraPlus™** skin barrier with **Remois Technology*** has a soft and flexible design to help minimize leakage with the added benefit of ceramide to help keep skin healthy. Ceramide is a natural component of human skin that helps prevent water loss that can lead to skin dryness and damage.



See pages
16 and 36
for more detail.



**Order a free
sample today.**

**Call your local Hollister
rep or visit**

www.hollister.com



Prior to use, be sure to read the Instructions for Use for information regarding Intended Use, Contraindications, Warnings, Precautions and Instructions for Use. The Hollister logo, CeraPlus and "Healthy skin. Positive outcomes." are trademarks of Hollister Incorporated. All other trademarks and copyrights are the property of their respective owners.
© 2020 Hollister Incorporated.



*Remois is a technology
of Alcare Co., Ltd.



Ostomy Care
Healthy skin. Positive outcomes.

New Image Two-Piece Pouching System

Drainable Pouches






1818x/1819x

Box of 10

- ▶ Drainable Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18192	18182	 1-3/4" (44 mm)	12" (30 cm)
18193	18183	 2-1/4" (57 mm)	12" (30 cm)
18194	18184	 2-3/4" (70 mm)	12" (30 cm)






1828x/1829x

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18292	18282	 1-3/4" (44 mm)	7" (18 cm)
18293	18283	 2-1/4" (57 mm)	7" (18 cm)
18294	18284	 2-3/4" (70 mm)	7" (18 cm)







1811x/1813x/1800x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

HCPCS: A5063

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18002	18132	18112	 1-3/4" (44 mm)	12" (30 cm)
18003	18133	18113	 2-1/4" (57 mm)	12" (30 cm)
18004	18134	18114	 2-3/4" (70 mm)	12" (30 cm)
18006	—	—	 4" (102 mm)	12" (30 cm)

Drainable Pouches (continued)






1814x/1816x

- ▶ Drainable Pouch
- ▶ With Integrated AF300 Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18162	18142	 1-3/4" (44 mm)	12" (30 cm)
18163	18143	 2-1/4" (57 mm)	12" (30 cm)
18164	18144	 2-3/4" (70 mm)	12" (30 cm)







1810x/1812x/1817x

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A5063

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18172	18102	18122	 1-3/4" (44 mm)	12" (30 cm)
18173	18103	18123	 2-1/4" (57 mm)	12" (30 cm)
18174	18104	18124	 2-3/4" (70 mm)	12" (30 cm)
18176	—	—	 4" (102 mm)	12" (30 cm)






1820x

- ▶ Drainable Mini-Pouch
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 10
Includes 1 clamp

HCPCS: A5063

Beige Stock No	Flange Size	Pouch Length
18202	 1-3/4" (44 mm)	9" (23 cm)
18203	 2-1/4" (57 mm)	9" (23 cm)
18204	 2-3/4" (70 mm)	9" (23 cm)



DID YOU KNOW...

- All *New Image* pouches are made up of multiple layers of odor-barrier film
- *New Image* barriers must be used with *New Image* pouches and should not be interchanged with other pouching systems

New Image Two-Piece Pouching System

High Output Drainable Pouches







1801x

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ Soft Drain Tap
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel — Body Side Only
- ▶ No adapter needed to connect to bedside collector

HCPCS: A4412

Ultra-Clear Stock No	Flange Size	Pouch Length
18012	 1-3/4" (44 mm)	12" (30 cm)
18013	 2-1/4" (57 mm)	12" (30 cm)
18014	 2-3/4" (70 mm)	12" (30 cm)
18016	 4" (102 mm)	12" (30 cm)

Closed Pouches






1836x/1837x

Box of 60

- ▶ Closed Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18362	18372	 1-3/4" (44 mm)	9" (23 cm)
18363	18373	 2-1/4" (57 mm)	9" (23 cm)
18364	18374	 2-3/4" (70 mm)	9" (23 cm)






1832x

Box of 30

- ▶ Closed Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *QuietWear* Pouch Material — Both Sides

HCPCS: A4419

Beige Stock No	Flange Size	Pouch Length
18322	 1-3/4" (44 mm)	9" (23 cm)
18323	 2-1/4" (57 mm)	9" (23 cm)
18324	 2-3/4" (70 mm)	9" (23 cm)

Closed Pouches (continued)



1838x/1839x

Box of 60

- ▶ Closed Mini-Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18382	18392	1-3/4" (44 mm)	7" (18 cm)
18383	18393	2-1/4" (57 mm)	7" (18 cm)
18384	18394	2-3/4" (70 mm)	7" (18 cm)



1833x

Box of 30

1873x

Box of 60

- ▶ Closed Pouch
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *QuietWear* Pouch Material — Both Sides

HCPCS: A5054

Beige Stock No Box of 30	Beige Stock No Box of 60	Flange Size	Pouch Length
18332	18732	1-3/4" (44 mm)	9" (23 cm)
18333	18733	2-1/4" (57 mm)	9" (23 cm)
18334	18734	2-3/4" (70 mm)	9" (23 cm)



1835x

Box of 30

1875x

Box of 60

- ▶ Closed Mini-Pouch
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A5054

Beige Stock No Box of 30	Beige Stock No Box of 60	Flange Size	Pouch Length
18352	18752	1-3/4" (44 mm)	7" (18 cm)
18353	18753	2-1/4" (57 mm)	7" (18 cm)
18354	—	2-3/4" (70 mm)	7" (18 cm)

Products for Colostomy Irrigations - see page 86-87

New Image Two-Piece Pouching System

Urostomy Pouches

Upgraded Design Recommended replacement for 1840x/1841x/1842xx series below



1890x/1891x/1892x




- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 2 adapters

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18922	18902	18912	 1-3/4" (44 mm)	9" (23 cm)
18923	18903	18913	 2-1/4" (57 mm)	9" (23 cm)
18924	18904	18914	 2-3/4" (70 mm)	9" (23 cm)

Original Design






1840x/1841x/1842x

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 adapter

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18422	18402	18412	 1-3/4" (44 mm)	9" (23 cm)
18423	18403	18413	 2-1/4" (57 mm)	9" (23 cm)
—	18404	18414	 2-3/4" (70 mm)	9" (23 cm)



DID YOU KNOW...

- The upgraded 1890x urostomy pouch series comes with two adapters. The original 1840x urostomy pouch series only comes with one
- Please see page 83 for the Hollister Bedside Drainage Collection System with Anti-Reflux Valve or visit www.hollister.com for all Hollister Leg Bag Systems

Non Sterile Kits — Drainable



1960x




Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear Pouch Film
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

Box of 5

Ceramide

HCPCS: A5063 + A4414

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19602	 1-3/4" (44 mm)	1-1/4" (32 mm)
19603	 2-1/4" (57 mm)	1-3/4" (44 mm)
19604	 2-3/4" (70 mm)	2-1/4" (57 mm)



1980x




Colostomy/Ileostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear Pouch Film
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only
- ▶ Clamp Included

Box of 5

Ceramide

HCPCS: A5063 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19802	 1-3/4" (44 mm)	1-1/4" (32 mm)
19803	 2-1/4" (57 mm)	1-3/4" (44 mm)
19804	 2-3/4" (70 mm)	2-1/4" (57 mm)



1900x





Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear Pouch Film
- ▶ (1) 7767 Loop Ostomy Bridge (in 19006 only)
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

Box of 5

HCPCS for 19002, 19003 & 19004: A5063 + A4414

HCPCS for 19006: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19002	 1-3/4" (44 mm)	1-1/4" (32 mm)
19003	 2-1/4" (57 mm)	1-3/4" (44 mm)
19004	 2-3/4" (70 mm)	2-1/4" (57 mm)
19006	 4" (102 mm)	3-1/2" (89 mm)



1910x





Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear Pouch Film
- ▶ (1) 7767 Loop Ostomy Bridge (in 19106 only)
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

Box of 5

HCPCS for 19102, 19103 & 19104: A5063 + A4414

HCPCS for 19106: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19102	 1-3/4" (44 mm)	1-1/4" (32 mm)
19103	 2-1/4" (57 mm)	1-3/4" (44 mm)
19104	 2-3/4" (70 mm)	2-1/4" (57 mm)
19106	 4" (102 mm)	3-1/2" (89 mm)

New Image Two-Piece Pouching System

Non Sterile Kits — Drainable (continued)





1930x

Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1410x *FormaFlex* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear
- ▶ (1) 1842x with *ComfortWear* Panels
- ▶ Single-Use Package

HCPCS: A5063 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19303	 2-1/4" (57 mm)	1-3/4" (44 mm)
19304	 2-3/4" (70 mm)	2-1/4" (57 mm)

Non Sterile Kits — Urostomy



1990x




Box of 5

Ceramide

Urostomy Kit includes:

- ▶ (1) 1120x *CeraPlus* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1892x with *ComfortWear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels-Body Side Only

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19902	 1-3/4" (44 mm)	1-1/4" (32 mm)
19903	 2-1/4" (57 mm)	1-3/4" (44 mm)
19904	 2-3/4" (70 mm)	2-1/4" (57 mm)






1920x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x *Flexend* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *Comfortwear* Panels-Body Side Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19202	 1-3/4" (44 mm)	1-1/4" (32 mm)
19203	 2-1/4" (57 mm)	1-3/4" (44 mm)
19204	 2-3/4" (70 mm)	2-1/4" (57 mm)





1940x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1410x *FormaFlex* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *Comfortwear* Panels-Body Side Only (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19402	 1-3/4" (44 mm)	1-1/4" (32 mm)
19403	 2-1/4" (57 mm)	1-11/16" (43 mm)

Sterile Kits — Drainable



1905x





Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear Pouch Film with *Comfortwear* Panels
- ▶ (1) 7767 Loop Ostomy Bridge (in 19056 only)
- ▶ Single-Use Package

HCPCS for 19052, 19053 & 19054: A5063 + A4414

HCPCS for 19056: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19052	 1-3/4" (44 mm)	1-1/4" (32 mm)
19053	 2-1/4" (57 mm)	1-3/4" (44 mm)
19054	 2-3/4" (70 mm)	2-1/4" (57 mm)
19056	 4" (102 mm)	3-1/2" (89 mm)



1915x





Box of 5
1 clamp per kit

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1420x *FlexWear* (standard wear) Skin Barrier with Tape Border
- ▶ (1) 1817x (Clamp Closure) Ultra-Clear with *Comfortwear* Panels
- ▶ (1) 7767 Loop Ostomy Bridge (in 19156 only)
- ▶ Single-Use Package

HCPCS for 19152, 19153 & 19154: A5063 + A4414

HCPCS for 19156: A5063 + A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19152	 1-3/4" (44 mm)	1-1/4" (32 mm)
19153	 2-1/4" (57 mm)	1-3/4" (44 mm)
19154	 2-3/4" (70 mm)	2-1/4" (57 mm)
19156	 4" (102 mm)	3-1/2" (89 mm)



1935x



Box of 5

Colostomy/Ileostomy Kit includes:

- ▶ (1) 1520x *FlexWear* (standard wear) Skin Barrier **without Tape Border**
- ▶ (1) 1800x (*Lock 'n Roll* Microseal Closure) Ultra-Clear with *Comfortwear* Panels
- ▶ Single-Use Package

HCPCS for 19353: A5063 & A4414

HCPCS for 19354: A5063 & A4415

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19353	 2-1/4" (57 mm)	1-3/4" (44 mm)
19354	 2-3/4" (70 mm)	2-1/4" (57 mm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

New Image Two-Piece Pouching System

Sterile Kits — Urostomy






1925x

Box of 5

Urostomy Kit includes:

- ▶ (1) 1460x *Flexend* (extended wear) Skin Barrier with Tape Border
- ▶ (1) 1842x with *ComfortWear* Panels (Ultra-Clear Pouch Film)
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 + A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19252	 1-3/4" (44 mm)	1-1/4" (32 mm)
19253	 2-1/4" (57 mm)	1-3/4" (44 mm)
19254	 2-3/4" (70 mm)	2-1/4" (57 mm)



1950x

Box of 5

Urostomy Kit includes:

- ▶ (1) 15603 *Flexend* (extended wear) Skin Barrier **without Tape Border**
- ▶ (1) 1842x with *ComfortWear* Panels
- ▶ (1) 18423 Ultra-Clear Pouch Film
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package

HCPCS: A4432 & A4409

Ultra-Clear Stock No	Flange Size	For Stomas Up To
19503	 2-1/4" (57 mm)	1-3/4" (44 mm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

A combination of the right

fit and formulation is the key

to positive outcomes



Premier
One-Piece Pouching System

Premier

One-Piece Pouching System

Why choose a one-piece pouching system?

With a one-piece pouching system, the skin barrier and pouch are connected and cannot be separated. Here are a few reasons why you might want to use a one-piece pouching system:

- You need a pouching system that is flexible
- You may prefer a pouching system that does not require putting pieces together
- You do not wish to change your pouch more often than your skin barrier

Upgraded Product Design

Select *Premier* One-Piece Pouching Systems feature two important product upgrades:



- **Lock 'n Roll** Microseal Closure with **ribbed strips** to facilitate ease of use, and **transparent flap** for increased security
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - **Flexible release liner** covers both the skin barrier and the tape border to help provide improved ease of use
 - **Tapered barrier edges** to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

Look for the **UPGRADED DESIGN** header above these stock numbers.

Go to page 57-59, 61-62 and 66-70 for more detail.

For more information, contact your Hollister Ostomy sales representative.



Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers



8901 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

Box of 10

Ceramide

HCPCS: A5056

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8901	2-1/2" (64 mm)	12" (30 cm)



8931 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 10

Ceramide

HCPCS: A4388

Ultra-Clear

Stock No	For Stomas Up To	Pouch Length
8931	2-1/2" (64 mm)	12" (30 cm)



8331/8371 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige with Viewing Option — Both Sides

Box of 10

HCPCS: A5056

Transparent Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8331	8371	2-1/2" (64 mm)	12" (30 cm)



8531/8551 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Beige — Both Sides
 - Ultra-Clear — Body Side Only

Box of 10

HCPCS: A4388

Ultra-Clear Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8531	8551	2-1/2" (64 mm)	12" (30 cm)

Drainable Pouches with Flat Barriers (continued)



8631/8651 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Includes 1 clamp

HCPCS: A4388

Ultra-Clear Stock No	Transparent Stock No	For Stomas Up To	Pouch Length
8651	8631	2-1/2" (64 mm)	12" (30 cm)



8641 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Mini-Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Includes 1 clamp

HCPCS: A4388

Transparent Stock No	For Stomas Up To	Pouch Length
8641	2-1/2" (64 mm)	9" (23 cm)



8081/8181/8381 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige and Beige with Viewing Option — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8081	8181	8381	2-1/2" (64 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



8031 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

8031	2-1/2" (64 mm)	12" (30 cm)
------	----------------	-------------



8131 (CUT-TO-FIT)

**Box of 10
Includes 1 clamp**

- ▶ Drainable Pouch
- ▶ *FlexWear* (standard wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

8131	2-1/2" (64 mm)	12" (30 cm)
------	----------------	-------------



88300/88400 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No

Beige Stock No

For Stomas Up To

Pouch Length

88400	88300	2-1/8" (55 mm)	12" (30 cm)
-------	-------	----------------	-------------



88800/88700 (CUT-TO-FIT)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No

Beige Stock No

For Stomas Up To

Pouch Length

88800	88700	2-1/8" (55 mm)	7" (18 cm)
-------	-------	----------------	------------

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Drainable Pouches with Flat Barriers (continued)



88302/88402 (OVAL CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ Without Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

OVAL

HCPCS: A4424

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
88402	88302	3" x 2.5" (75 mm x 65 mm)	12" (30 cm)



88401 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
88401	2-1/2" (64 mm)	12" (30 cm)



88501 (CUT-TO-FIT)

Box of 10
Includes 1 clamp

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ Without Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

Ultra-Clear Stock No	For Stomas Up To	Pouch Length
88501	2-1/2" (64 mm)	12" (30 cm)



89xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

Ceramide

HCPCS: A5056

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8925	1" (25 mm)	12" (30 cm)
8930	1-3/16" (30 mm)	12" (30 cm)
8935	1-3/8" (35 mm)	12" (30 cm)
8940	1-9/16" (40 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Flat Barriers (continued)



853x/855x (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4388

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
—	8552	1" (25 mm)	12" (30 cm)
8538	8558	1-1/4" (32 mm)	12" (30 cm)



863x (PRE-SIZED)

Box of 10
Includes 1 clamp

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

Transparent Stock No	Barrier Opening	Pouch Length
8637	3/4" (19 mm)	12" (30 cm)
8632	1" (25 mm)	12" (30 cm)
8638	1-1/4" (32 mm)	12" (30 cm)
8633	1-1/2" (38 mm)	12" (30 cm)
8639	1-3/4" (44 mm)	12" (30 cm)
8634	2" (51 mm)	12" (30 cm)
8635	2-1/2" (64 mm)	12" (30 cm)



864x (PRE-SIZED)

Box of 10
Includes 1 clamp

- ▶ Drainable Mini-Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A4388

Beige Stock No	Barrier Opening	Pouch Length
8647	3/4" (19 mm)	9" (23 cm)
8642	1" (25 mm)	9" (23 cm)
8648	1-1/4" (32 mm)	9" (23 cm)
8643	1-1/2" (38 mm)	9" (23 cm)
8649	1-3/4" (44 mm)	9" (23 cm)

Drainable Pouches with Flat Barriers (continued)



883xx/884xx (PRE-SIZED)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4424

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
88425	88325	1" (25 mm)	12" (30 cm)
88430	88330	1-3/16" (30 mm)	12" (30 cm)
88435	88335	1-3/8" (35 mm)	12" (30 cm)
88440	88340	1-9/16" (40 mm)	12" (30 cm)



887xx (PRE-SIZED)

Box of 20

- ▶ Drainable Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A4424

Beige Stock No	Barrier Opening	Pouch Length
88725	1" (25 mm)	7" (18 cm)
88730	1-3/16" (30 mm)	7" (18 cm)
88735	1-3/8" (35 mm)	7" (18 cm)
88740	1-9/16" (40 mm)	7" (18 cm)



A flexible fit


made secure

Premier Soft Convex Pouching System

Designed to stay secure, no matter what shape the body takes.

When it comes to finding the right solution for patients, fit comes first. Comfort and ease-of-use are also essential to help ensure positive outcomes—in the hospital and after patients leave your care.

Hollister Ostomy presents Premier soft convex—a pouching system that delivers both convex support and flexibility. The skin barrier is designed to conform to the body, without compromising fit. Belt tabs and a tape-bordered barrier help provide an extra sense of security. And the new viewing option offers a convenient way to help position the pouch and observe output.

- Look for the soft convex symbol 
- Available in CeraPlus and Flextend (extended wear) skin barrier formulations

Premier One-Piece Pouching System

Drainable Pouches with Soft Convex Barriers



895xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option

Stock No	For Stomas Up To	Pouch Length
8954	1" (25 mm)	12" (30 cm)
8958	1-1/2" (38 mm)	12" (30 cm)
89511	2-1/8" (55 mm)	12" (30 cm)



857xx/867xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *Flextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige with Viewing Option — Both Sides

Box of 5



HCPCS: A5057

Ultra-Clear Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8574	8674	1" (13-25 mm)	12" (30 cm)
8578	8678	1 1/2" (13-38 mm)	12" (30 cm)
85711	86711	2 1/8" (13-55 mm)	12" (30 cm)



896x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option Stock No	Opening	Pouch Length
8960	3/4" (19 mm)	12" (30 cm)
8961	7/8" (22 mm)	12" (30 cm)
8962	1" (25 mm)	12" (30 cm)
8963	1-1/8" (29 mm)	12" (30 cm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Premier One-Piece Pouching System

Drainable Pouches with Soft Convex Barriers (continued)



866x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ Soft Convex
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option Stock No

Barrier Opening

Pouch Length

8661	7/8" (22 mm)	12" (30 cm)
8662	1" (25 mm)	12" (30 cm)
8663	1-1/8" (29 mm)	12" (30 cm)
8664	1-1/4" (32 mm)	12" (30 cm)
8665	1-1/2" (38 mm)	12" (30 cm)
8666	1-3/4" (44 mm)	12" (30 cm)

Drainable Pouches with Convex Barriers



891xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS Code: A5057

Beige w/Viewing Option Stock No

For Stomas Up To

Pouch Length

8914	1" (25 mm)	12" (30 cm)
8918	1-1/2" (38 mm)	12" (30 cm)
89111	2" (51 mm)	12" (30 cm)



898xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5



HCPCS Code: A4390

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

8984	1" (25 mm)	12" (30 cm)
8988	1-1/2" (38 mm)	12" (30 cm)
89811	2" (51 mm)	12" (30 cm)

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 858xx series below



827xx/854x (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige with Viewing Option — Both Sides

Box of 5



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS Code: A5057

Transparent Stock No	Beige w/Viewing Option Stock No	For Stomas Up To	Pouch Length
8274	8544	1" (25 mm)	12" (30 cm)
8278	8548	1-1/2" (38 mm)	12" (30 cm)
82711	—	2" (51 mm)	12" (30 cm)

Original Design

858xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only
- ▶ *Lock 'n Roll* Closure

Box of 5



HCPCS Code: A5057

Transparent Stock No	For Stomas Up To	Pouch Length
8584	1" (25 mm)	12" (30 cm)
8588	1-1/2" (38 mm)	12" (30 cm)
85811	2" (51 mm)	12" (30 cm)



DID YOU KNOW...

Select **Premier** one-piece pouching systems feature two important product upgrades:

- *Lock 'n Roll* Microseal Closure with ribbed strips to facilitate ease of use, and transparent flap for increased security
- Upgraded **skin barrier design** to help increase comfort and help improve ease of use:
 - Flexible release liner covers both the skin barrier and the tape border to help provide improved ease of use
 - Tapered barrier edges to help provide increased comfort
 - No change in skin barrier application
 - Same HCPCS code

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 850xx/852xx series below



824xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only

Box of 5



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HPCPS: A4390

Transparent Stock No	For Stomas Up To	Pouch Length
8244	1" (25 mm)	12" (30 cm)
8248	1-1/2" (38 mm)	12" (30 cm)
82411	2" (51 mm)	12" (30 cm)

Original Design

850xx/852xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides
- ▶ *Lock 'n Roll* Closure

Box of 5



HPCPS: A4390

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8524	8504	1" (25 mm)	12" (30 cm)
8528	8508	1-1/2" (38 mm)	12" (30 cm)
85211	85011	2" (51 mm)	12" (30 cm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 862xx series below



862xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
Includes 1 clamp



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HPCPS: A4390

Transparent Stock No

For Stomas Up To

Pouch Length

82624	1" (25 mm)	12" (30 cm)
82628	1-1/2" (38 mm)	12" (30 cm)
826211	2" (51 mm)	12" (30 cm)

Original Design

862xx (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
Includes 1 clamp



HPCPS: A4390

Transparent Stock No

For Stomas Up To

Pouch Length

8624	1" (25 mm)	12" (30 cm)
8628	1-1/2" (38 mm)	12" (30 cm)
86211	2" (51 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)



899x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Both Sides

Box of 5



HCPCS: A5057

Beige w/Viewing Option

Stock No	Barrier Opening	Pouch Length
8990	3/4" (19 mm)	12" (30 cm)
8991	7/8" (22 mm)	12" (30 cm)
8992	1" (25 mm)	12" (30 cm)
8993	1-1/8" (29 mm)	12" (30 cm)



DID YOU KNOW...

- *Premier* cut-to-fit products have a 1/2" (13 mm) starter hole and can be used as a pre-sized product for some small stomas (excludes high output pouches, wound drainage collectors and *Pouchkins* products)
- Cutting guides are provided for all cut-to-fit barriers
- Drainable pouches without an integrated closure include one pouch clamp per box
- We do not recommend protective wipes be used with extended wear skin barriers, as they may decrease wear time
- Filtered pouches do not require a cover or sticker when the pouch may get wet
- All convex *Premier* pouches have tape border barriers
- All convex *Premier* pouches have belt tabs

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 851xx/859x series below



828x/829x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides



Upgraded Design Features:

- ▶ *Lock 'n Roll* Microseal Closure with Ribbed Strips and Transparent Flap
- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8280	8290	3/4" (19 mm)	12" (30 cm)
8281	8291	7/8" (22 mm)	12" (30 cm)
8282	8292	1" (25 mm)	12" (30 cm)
8283	—	1-1/8" (29 mm)	12" (30 cm)
8284	8294	1-1/4" (32 mm)	12" (30 cm)
8285	8295	1-3/8" (35 mm)	12" (30 cm)
8286	—	1-1/2" (38 mm)	12" (30 cm)
8288	—	1-3/4" (44 mm)	12" (30 cm)

Original Design

851xx/859x (PRE-SIZED)

Box of 5

- ▶ Drainable Pouch
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides
- ▶ *Lock 'n Roll* Closure



HCPCS: A4390

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8510	8590	3/4" (19 mm)	12" (30 cm)
8511	8591	7/8" (22 mm)	12" (30 cm)
8512	8592	1" (25 mm)	12" (30 cm)
8513	8593	1-1/8" (29 mm)	12" (30 cm)
8514	8594	1-1/4" (32 mm)	12" (30 cm)
8515	8595	1-3/8" (35 mm)	12" (30 cm)
85116	8596	1-1/2" (38 mm)	12" (30 cm)
—	8597	1-5/8" (41 mm)	12" (30 cm)
85118	—	1-3/4" (44 mm)	12" (30 cm)
85119	—	2" (51 mm)	12" (30 cm)

Premier One-Piece Pouching System

Drainable Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 861xx series below



8261x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
Includes 1 clamp



Upgraded Design Features:

- ▶ Flexible Release Liner
- ▶ Tapered Barrier Edges

HCPCS: A4390

Transparent Stock No

Transparent Stock No	Barrier Opening	Pouch Length
82610	3/4" (19 mm)	12" (30 cm)
82611	7/8" (22 mm)	12" (30 cm)
82612	1" (25 mm)	12" (30 cm)
82613	1-1/8" (29 mm)	12" (30 cm)
82614	1-1/4" (32 mm)	12" (30 cm)
82615	1-3/8" (35 mm)	12" (30 cm)
82616	1-1/2" (38 mm)	12" (30 cm)
82617	1-5/8" (41 mm)	12" (30 cm)
82619	2" (51 mm)	12" (30 cm)

Original Design

861x (PRE-SIZED)

- ▶ Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 5
Includes 1 clamp



HCPCS: A4390

Transparent Stock No

Transparent Stock No	Barrier Opening	Pouch Length
8610	3/4" (19 mm)	12" (30 cm)
8611	7/8" (22 mm)	12" (30 cm)
8612	1" (25 mm)	12" (30 cm)
8613	1-1/8" (29 mm)	12" (30 cm)
8614	1-1/4" (32 mm)	12" (30 cm)
8615	1-3/8" (35 mm)	12" (30 cm)
8616	1-1/2" (38 mm)	12" (30 cm)
8617	1-5/8" (41 mm)	12" (30 cm)
8618	1-3/4" (44 mm)	12" (30 cm)
8619	2" (51 mm)	12" (30 cm)

High Output Drainable Pouches



80070/80110 (CUT-TO-FIT)

Box of 10

- ▶ Non Sterile, Drainable Pouch
- ▶ *Flexextend* (extended wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Soft Tap
- ▶ With Access Window
- ▶ With *ComfortWear* Panels — Body Side Only
- ▶ Without Starter Hole
- ▶ No adapter needed to connect to bedside collector

HCPCS: A4435

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

80070	2-3/4" (70 mm)	12" (30 cm)
80110	4-1/3" (110 mm)	12" (30 cm)

Wound Drainage Collector



970x/977x

- ▶ For Management of Draining Wounds, Fistulae, or Draining Tube Sites
- ▶ Available With and Without *FlexWear* (standard wear) Skin Barriers — Flat
- ▶ **Without Tape Border**
- ▶ Sterile and Non Sterile Drainable Pouches
- ▶ With Access Window
- ▶ Reclosable Cap
- ▶ Without Starter Hole
- ▶ No adapter needed to connect to bedside collector

HCPCS: A6154

With Barrier: Non Sterile

Stock No	Description	Quantity/Box
9773	Small (for wounds up to 3")	10
9776	Medium (for wounds up to 3-3/4")	10
9778	Large (for wounds up to 4" x 8")	3

With Barrier: Sterile

Stock No	Description	Quantity/Box
9701	Small (for wounds up to 3")	5
9702	Medium (for wounds up to 3-3/4")	5
9703	Large (for wounds up to 4" x 8")	3

Without Barrier: Sterile

Stock No	Description	Quantity/Box
9772	Small (for wounds up to 3")	10
9775	Medium (for wounds up to 3-3/4")	10



DID YOU KNOW...

- All high output pouches and wound drainage collectors can be connected directly into a bedside collector without an adaptor

Premier One-Piece Pouching System

Closed Pouches with Flat Barriers



82300/82400 (CUT-TO-FIT)

Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82400	82300	2-1/8" (55 mm)	9" (23 cm)



82302/82402 (OVAL CUT-TO-FIT)

Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

OVAL

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82402	82302	3" x 2.5" (75 mm x 65 mm)	9" (23 cm)



82100/82500 (CUT-TO-FIT)

Box of 30

- ▶ Closed Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4416

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
82500	82100	2-1/8" (55 mm)	7" (18 cm)



823xx (PRE-SIZED)

Box of 30

- ▶ Closed Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels — Both Sides

HCPCS: A4416

Beige Stock No	Barrier Opening	Pouch Length
82325	1" (25 mm)	9" (23 cm)
82330	1-3/16" (30 mm)	9" (23 cm)
82335	1-3/8" (35 mm)	9" (23 cm)

Closed Pouches with Flat Barriers (continued)



821xx (PRE-SIZED)

Box of 30

- ▶ Closed Mini-Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated *AF300* Filter
- ▶ With *ComfortWear* Panels — Both Sides

HC416: A4416

Beige Stock No	Barrier Opening	Pouch Length
82125	1" (25 mm)	7" (18 cm)
82130	1-3/16" (30 mm)	7" (18 cm)
82135	1-3/8" (35 mm)	7" (18 cm)

Stoma Cap



1796 (PRE-SIZED)

Box of 30

- ▶ Stoma Cap
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ With Integrated Filter

HCPCS: A5055

Beige Stock No	Barrier Opening	Size
1796	1-15/16" (50 mm)	4" (102 mm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Premier One-Piece Pouching System

Urostomy Pouches with Flat Barriers

Upgraded Design Recommended replacement for 8440/8450/8460 series below



84490/84590/84690 (CUT-TO-FIT) **Box of 10**
Includes 2 adapters

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Tapered Barrier Edges
- ▶ Drain Valve Indicator

HCPCS: A4428

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
84590	84690	84490	2-1/2" (64 mm)	9" (23 cm)

Original Design



8440/8450/8460 (CUT-TO-FIT) **Box of 10**
Includes 1 adapter

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra Clear — Body Side Only
 - Beige — Both Sides

HCPCS: A4428

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8450	8460	8440	2-1/2" (64 mm)	9" (23 cm)

The upgraded Hollister urostomy pouch is available in the **Premier** one-piece pouching system.

The upgraded pouch features several innovations designed to deliver more comfort and a greater sense of control:

- Skin barrier options to create custom solutions for your patients' unique needs
 - ▶ Flat, convex or soft convex
 - ▶ Tape border
 - ▶ *Flexend* (extended wear) skin barrier
 - ▶ *CeraPlus* (extended wear) skin barrier
- Enhanced flow-control tap is made with soft material and a streamlined tap
- Multi-chamber design helps to facilitate even distribution of urine
- Improved adapter helps facilitate ease of use
- Belt tab option (convex)

Urostomy Pouches with Flat Barriers (continued)

Upgraded Design Recommended replacement for 844x/846x series below



8449x/8469x (PRE-SIZED)

Box of 10
Includes 2 adapters

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Tapered Barrier Edges
- ▶ Drain Valve Indicator

HCPCS: A4428

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
—	84497	3/4" (19 mm)	9" (23 cm)
84692	—	1" (25 mm)	9" (23 cm)
84698	84498	1-1/4" (32 mm)	9" (23 cm)
84693	—	1-1/2" (38 mm)	9" (23 cm)
84699	84499	1-3/4" (44 mm)	9" (23 cm)
84694	84494	2" (51 mm)	9" (23 cm)

Original Design



844x/846x (PRE-SIZED)

Box of 10
Includes 1 adapter

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4428

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8467	8447	3/4" (19 mm)	9" (23 cm)
8462	8442	1" (25 mm)	9" (23 cm)
8468	8448	1-1/4" (32 mm)	9" (23 cm)
8463	8443	1-1/2" (38 mm)	9" (23 cm)
8469	8449	1-3/4" (44 mm)	9" (23 cm)
8464	—	2" (51 mm)	9" (23 cm)

Premier One-Piece Pouching System

Urostomy Pouches with Convex Barriers

Upgraded Design Recommended replacement for 847x/843x series below



8479xxx/8439xxx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 5
Includes 1 adapter



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Tapered Barrier Edges
- ▶ Drain Valve Indicator

HCPCS: A4430

Ultra-Clear Stock No	Beige Stock No	For Stomas Up To	Pouch Length
84794	—	1" (25 mm)	9" (23 cm)
84798	—	1-1/2" (38 mm)	9" (23 cm)
8479111	8439111	2" (51 mm)	9" (23 cm)

Original Design



847xx/843xx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 5
Includes 1 adapter



HCPCS: A4430

Transparent Stock No	Beige Stock No	For Stomas Up To	Pouch Length
8474	—	1" (25 mm)	9" (23 cm)
8478	—	1-1/2" (38 mm)	9" (23 cm)
84711	—	2" (51 mm)	9" (23 cm)



DID YOU KNOW...

- Leaving urine in the urostomy pouch before connecting to a bedside collector/leg bag may minimize the suction effect and allow urine to drain better into the collection system

Urostomy Pouches with Convex Barriers (continued)

Upgraded Design Recommended replacement for 848x/849x series below



8489x/8499x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 5
Includes 1 adapter



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Tapered Barrier Edges
- ▶ Drain Valve Indicator

HCPCS: A4430

Ultra Clear Stock No	Beige Stock No	Barrier Opening	Pouch Length
84890	—	1/2" (13 mm)	9" (23 cm)
84891	84991	5/8" (16 mm)	9" (23 cm)
84892	84992	3/4" (19 mm)	9" (23 cm)
84893	84993	7/8" (22 mm)	9" (23 cm)
84894	84994	1" (25 mm)	9" (23 cm)
84895	84995	1-1/8" (29 mm)	9" (23 cm)
84896	84996	1-1/4" (32 mm)	9" (23 cm)
84897	84997	1-3/8" (35 mm)	9" (23 cm)
84898	84998	1-1/2" (38 mm)	9" (23 cm)
84899	84999	1-3/4" (44 mm)	9" (23 cm)

Original Design



848x/849x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

Box of 5
Includes 1 adapter



HCPCS: A4430

Transparent Stock No	Beige Stock No	Barrier Opening	Pouch Length
8480	—	1/2" (13 mm)	9" (23 cm)
8481	8491	5/8" (16 mm)	9" (23 cm)
8482	8492	3/4" (19 mm)	9" (23 cm)
8483	8493	7/8" (22 mm)	9" (23 cm)
8484	8494	1" (25 mm)	9" (23 cm)
8485	8495	1-1/8" (29 mm)	9" (23 cm)
8486	8496	1-1/4" (32 mm)	9" (23 cm)
8487	8497	1-3/8" (35 mm)	9" (23 cm)
8488	8498	1-1/2" (38 mm)	9" (23 cm)
8489	—	1-3/4" (44 mm)	9" (23 cm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Premier One-Piece Pouching System

Urostomy Pouches with **CeraPlus** Soft Convex Barriers

Upgraded Design



8413xx (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
includes 1 adapter**



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator
- ▶ Tapered barrier edges

HCPCS: A4430

Ultra Clear Stock No	For Stomas Up To	Pouch Length
84134	1" (25 mm)	9" (23 cm)
84138	1-1/2" (38 mm)	9" (23 cm)
841311	2-1/8" (55 mm)	9" (23 cm)

Upgraded Design



841x (PRE-SIZED)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *CeraPlus* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels — Body Side Only

**Box of 5
includes 1 adapter**



Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator
- ▶ Tapered barrier edges

HCPCS: A4430

Beige w/Viewing Option Stock No	Barrier Opening	Pouch Length
8412	3/4" (19 mm)	9" (23 cm)
8413	7/8" (22 mm)	9" (23 cm)
8414	1" (25 mm)	9" (23 cm)
8415	1-1/8" (29 mm)	9" (23 cm)

Non Sterile Kits — Drainable



89001 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88401 *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89001

2-1/2" (64 mm)

12" (30 cm)



89002 (CUT-TO-FIT)

**Box of 5 Kits
1 clamp per kit**

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 88501 *SoftFlex* (standard wear) Skin Barrier
- ▶ **Without Tape Border**
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A5061

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89002

2-1/2" (64 mm)

12" (30 cm)



89003 (CUT-TO-FIT)

**Box of 5 Kits
1 clamp per kit**

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 8651 *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ Clamp Closure
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89003

2-1/2" (64 mm)

12" (30 cm)



89004 (CUT-TO-FIT)

Box of 5 Kits

Colostomy/Ileostomy Kit Includes:

- ▶ (1) 8531 *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4388

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89004

2-1/2" (64 mm)

12" (30 cm)

Non Sterile Kits — Urostomy



89005 (CUT-TO-FIT)

**Box of 5 Kits
1 adapter per kit**

Urostomy Kit Includes:

- ▶ (1) 8450 *Flexextend* (extended wear) Skin Barrier
- ▶ With Tape Border
- ▶ (1) 7331 Urostomy Drain Tube Adapter
- ▶ Single-Use Package
- ▶ With *ComfortWear* Panels — Body Side Only

HCPCS: A4428

**Ultra-Clear
Stock No**

For Stomas Up To

Pouch Length

89005

2-1/2" (64 mm)

9" (23 cm)



A secure fit is necessary for

preventing leakage and

providing patient comfort

Pouchkins
Pediatric Pouching System

Pouchkins
Pediatric Pouching System

Pouchkins
Pediatric Pouching System

Pouchkins Pediatric Pouching System

Premie Pouch



3777 (CUT-TO-FIT)

Box of 15

- ▶ Use as a Premie Closed or Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ Without Starter Hole
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure (if drain opening is desired)

Fun Character Stickers Included!

HCPCS: A5052 + A4362

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

3777

5/8" (15 mm)

4-1/2" (11.5 cm)

Newborn Pouch



3778 (CUT-TO-FIT)

Box of 15

- ▶ Newborn Pouch; Use With or Without the Drain Valve
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ Without Starter Hole
- ▶ Without Integrated Filter
- ▶ Soft Wire Closure (if larger drain opening is desired)

Fun Character Stickers Included!

HCPCS: A5062 + A4362

Ultra-Clear Stock No

For Stomas Up To

Pouch Length

3778

1-3/8" x 7/8" (35 mm x 22 mm)

6-1/2" (16 cm)

One-Piece Pouches



3796 (CUT-TO-FIT)

Box of 10

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ **Without Tape Border**
- ▶ Without Starter Hole
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Fun Character Stickers Included!

HCPCS: A5061

Transparent Stock No

For Stomas Up To

Pouch Length

3796

1-1/2" (38 mm)

6-1/2" (16 cm)

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Pouchkins Pediatric Pouching System

One-Piece Pouches (continued)



3795 (CUT-TO-FIT)

- ▶ Drainable Pouch
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ With Integrated AF300 Filter
- ▶ Lock 'n Roll Microseal Closure
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 10

Fun Character Stickers Included!

HCPCS: A4424

Transparent Stock No

For Stomas Up To

Pouch Length

3795

2" (51 mm)

7" (18 cm)



3797 (CUT-TO-FIT)

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border
- ▶ With *ComfortWear* Panels — Body Side Only

Box of 10

Includes 1 adapter

Fun Character Stickers Included!

HCPCS: A5071

Transparent Stock No

For Stomas Up To

Pouch Length

3797

1-1/2" (38 mm)

8-3/4" (22 cm)



DID YOU KNOW...

Pouchkins *SoftFlex* skin barrier

- Designed to be gentle
- No starter hole, which enables the opening to be placed where needed within the cutting guide
- No tape by design
- Flexible and sized for babies and children
- None of the *Pouchkins* pediatric pouching systems are made with natural rubber latex
- Any 1-3/4" *New Image* pouch will fit this skin barrier



Two-Piece Skin Barrier




3761 (CUT-TO-FIT)

Box of 5

- ▶ *SoftFlex* (standard wear) Skin Barrier — Flat
- ▶ Without Starter Hole
- ▶ Without Tape Border

HCPCS: A4414

Transparent Stock No	Flange Size	For Stomas Up To
3761	 1-3/4" (44 mm)	1-1/4" (32 mm)

Two-Piece Pouches



3799

Box of 10

- ▶ Drainable Pouch
- ▶ Without Integrated Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel — Body Side Only

HCPCS: A5063

Transparent Stock No	Flange Size	Pouch Length
3799	 1-3/4" (44 mm)	6-1/2" (16 cm)




18292/18282

Box of 20

- ▶ Drainable Pouch
- ▶ With Integrated *AF300* Filter
- ▶ *Lock 'n Roll* Microseal Closure
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panel
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4425

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18292	18282	 1-3/4" (44 mm)	7" (18 cm)



18382/18392

Box of 60

- ▶ Closed Mini-Pouch
- ▶ With Integrated *AF300* Filter
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent — Body Side Only
 - Beige — Both Sides

HCPCS: A4419

Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18382	18392	 1-3/4" (44 mm)	7" (18 cm)

Pouchkins Pediatric Pouching System

Two-Piece Pouches (continued)

Upgraded Design Recommended replacement for 1840x/1841x/1842xx series below



1890x/1891x/1892x

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 2 adapters

Upgraded Design Features:

- ▶ Multi-Chamber Design
- ▶ Softer, Streamlined Tap
- ▶ Drain Valve Indicator

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18922	18902	18912	 1-3/4" (44 mm)	9" (23 cm)

Original Design



18402/18412/18422

- ▶ Urostomy Pouch
- ▶ Anti-Reflux Valve
- ▶ With Belt Tabs
- ▶ With *ComfortWear* Panels
 - Transparent and Ultra-Clear — Body Side Only
 - Beige — Both Sides

Box of 10
Includes 1 adapter

HCPCS: A4432

Ultra-Clear Stock No	Transparent Stock No	Beige Stock No	Flange Size	Pouch Length
18422	18402	18412	 1-3/4" (44 mm)	9" (23 cm)

Pediatric Ostomy Belt



3774

- ▶ Ostomy Belt
- White Color

Single Unit

HCPCS: A4367

Stock No	Length
3774	10" - 17" (25 cm-43 cm) Adjustable



DID YOU KNOW...

Shadow Buddies dolls

- Exclusive to **Hollister**, *Shadow Buddies* dolls are special friends for children who have had or will have ostomy surgery
- Great teaching tool to help with education and care giving
- Huggable dolls with heart-shaped eyes and a stoma provide comfort and support



To order, see page 93.

Shadow Buddies is a trademark of The Shadow Buddies Foundation

Accessories can enhance

performance of the skin barrier

and promote healthy skin



Accessories

Accessories

Why use an ostomy accessory?

Ostomy accessories are used together with your pouching system. You may use one or more of these products. Here are a few reasons why you might want to use an ostomy accessory:

- It may help enhance the performance of your pouching system. For example, an accessory could help increase wear time, protect your skin, or make the removal of your adhesive easier
- It may help solve a problem such as leakage or skin irritation
- It may help improve the quality of your life by eliminating odor

Count on the Right Fit With

Adapt Barrier Rings



Hollister understands that no two people or stomas are alike — and that barrier rings can play a critical role in helping ensure the best outcome for your patients

Whether you choose flat or convex, slim or regular, round or oval—
Adapt barrier rings:

- Are designed to be stretched or shaped or cut into pieces to use as a filler in uneven areas
- Help prevent stoma output from getting under the pouch seal, which may reduce leakage and skin irritation
- Can be cut or stacked together to improve the fit of the pouching system
- Are alcohol free, so there is no sting from alcohol when applied to broken or irritated skin

Talk to your Hollister ostomy representative about the full line of Adapt barrier rings, now available with and without ceramide.



Ostomy Care
Healthy skin. Positive outcomes.

Adapt Skin Barrier Rings



88xx

- ▶ CeraRing Barrier Rings (Flat)
- ▶ Flextend M (extended wear) with Ceramide

Box of 10

Ceramide

HCPCS: A4385

Stock No	Outer Diameter	Thickness
8805	2" (48 mm)	4.5 mm
8815	2" (48 mm)	2.3 mm SLIM



78xx

- ▶ Barrier Rings (Flat)
- ▶ Flextend M (extended wear)

Box of 10

HCPCS: A4385

Stock No	Outer Diameter	Thickness
7805	2" (48 mm)	4.5 mm
7815	2" (48 mm)	2.3 mm SLIM
7806	4" (98 mm)	2.3 mm SLIM



895xx

- ▶ CeraRing Round Convex Barrier Rings
- ▶ Flextend (extended wear) with Ceramide

Box of 10

CONVEX
Ceramide

HCPCS: A4411

Stock No	Inner Diameter
89520	13/16" (20 mm) — can be stretched to 1" (25 mm)
89530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
89540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)



795xx

- ▶ Convex Round Convex Barrier Rings
- ▶ Flextend (extended wear)

Box of 10

CONVEX

HCPCS: A4411

Stock No	Inner Diameter
79520	13/16" (20 mm) — can be stretched to 1" (25 mm)
79530	1-3/16" (30 mm) — can be stretched to 1-3/8" (35 mm)
79540	1-9/16" (40 mm) — can be stretched to 1-3/4" (45 mm)



896xx

- ▶ CeraRing Oval Convex Barrier Rings
- ▶ Flextend (extended wear) with Ceramide

Box of 10

CONVEX
Ceramide

HCPCS: A4411

Stock No	Inner Diameter
89601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
89602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
89603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)



796xx

- ▶ Oval Convex Barrier Rings
- ▶ Flextend (extended wear)

Box of 10

CONVEX

HCPCS: A4411

Stock No	Inner Diameter
79601	7/8" x 1-1/2" (22 mm x 38 mm) — can be stretched to 1-1/8" x 1-3/4" (27 mm x 43 mm)
79602	1-3/16" x 1-7/8" (30 mm x 48 mm) — can be stretched to 1-3/8" x 2-1/8" (35 mm x 53 mm)
79603	1-1/2" x 2-3/16" (38 mm x 56 mm) — can be stretched to 1-3/4" x 2-3/8" (43 mm x 61 mm)



A Sense of Security

Means Everything



Adapt Barrier Extenders

At Hollister Ostomy Care, we know a sense of security and leakage prevention are top priorities for you.

Three strips create a custom frame around your skin barrier.

Our Adapt barrier extenders offer an additional layer of protection for those concerned about lifting around the edges of their skin barrier. With a thin and flexible hydrocolloid formulation, Adapt barrier extenders conform easily to body contours. Three convenient strips make it easy to create a custom frame around the skin barrier. The end result is enhanced security, when it matters most.

Ostomy Accessories

Adapt Skin Barrier Strips



79400

- ▶ Skin Barrier Strips
- ▶ *Flexextend* (extended wear)

Box of 10

HCPCS: A4385

Stock No	Size	Length
79400	60 g	5-1/8" (13 cm)

Adapt Barrier Extenders



79402

- ▶ Barrier Extenders
- ▶ Hydrocolloid Formulation

Box of 30

HCPCS: A4362*

Stock No
79402

* 3 extenders = 1 unit, Box of 30 = 10 billable units

Adapt Skin Barrier Paste



79300/79301

- ▶ Non-Adhesive

HCPCS: A4406

Stock No	Description	Quantity
79300	2.1 oz (60 g) tube	Box of 1
79301	0.5 oz (14 g) tube	Box of 20

Adapt Stoma Powder



7906

Box of 1

HCPCS: A4371

Stock No	Description
7906	1 oz (28.3 g) puff bottle

Adapt Odor Eliminator and Lubricant



78500/78501

HCPCS Code: A4394

Stock No	Description	Quantity
78500	8 oz (236 mL) bottle	Box of 1
78501	0.27 oz (8 mL) packet	50 Packets/Box

Adapt Ostomy Belts



7300/7299

Box of 10

HCPCS: A4367

Beige Stock No	Size
7300	23"-43" (58 cm-109 cm) Medium - Adjustable
7299	34"-65" (86 cm-165 cm) Large - Adjustable

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Adapt Skin Prepping Product



7917

- ▶ Skin Protective Wipes
- ▶ Non-Adhesive
- ▶ Alcohol Free

**1 Box
50 Wipes Each**

NO STING

HCPCS: A5120

Stock No

7917

**Not recommended for use with extended wear skin barriers, as it may decrease wear time*

Adapt Adhesives and Removers



7760

- ▶ Universal Remover Wipes
- ▶ For Adhesives and Barriers
- ▶ Silicone Based
- ▶ Alcohol Free

**1 Box
50 Wipes Each**

NO STING

HCPCS: A4456

Stock No

7760



7737

- ▶ Medical Adhesive Remover
- ▶ Silicone Based
- ▶ Alcohol Free

Box of 1

NO STING

HCPCS: A4455

Stock No

Description

7737

1.7 oz (50 mL) 360° spray can



7730

- ▶ Medical-Grade, Silicone-Based Adhesive
- ▶ Improves Adhesive Contact Between Skin and Barrier

Box of 1

HCPCS: A4364

Stock No

Description

7730

3.8 oz (112 mL) 360° spray can

Flexlend (extended wear) Skin Barrier



8800

8801

HCPCS for 8800: A4385

HCPCS for 8801: A5122

Stock No

Size

8800

4 x 4" (10 cm x 10 cm)

Box of 5

8801

8 x 8" (20 cm x 20 cm)

Box of 3

Ostomy Accessories

Odor Eliminators



773x

- ▶ *m9* Odor Eliminator Spray
- ▶ Neutralizes Odor Rather Than Masking It

HCPCS: No Code

Stock No (unscented)	Description	Quantity
7732	2 oz (60 mL) pump spray	Box of 12
7733	8 oz (240 mL) pump spray	Box of 6

(Green Apple Scented)

Stock No	Description	Quantity
7734	2 oz (60 mL) pump spray	Box of 12
7735	8 oz (240 mL) pump spray	Box of 6



771x

- ▶ *m9* Odor Eliminator Drops

HCPCS: A4394

Stock No	Description	Quantity
7715	1 oz (30 mL) bottle	Box of 12
7717	8 oz (240 mL) bottle	Box of 6

Cleaner/Decrystallizer



7736

- ▶ *m9* Cleaner/Decrystallizer
- ▶ Cleans Urinary Drainage Systems
- ▶ pH Balanced and Non-Acidic
- ▶ Free Wash Bottle Included

Box of 12

HCPCS: A5131

Stock No	Description
7736	16 oz (480 mL) bottle

Drainable Pouch Clamp



8770

HCPCS: A4363

Box of 20

Stock No	Color
8770	Beige

Urostomy Drain Tube Adapter



7331

HCPCS: A4421

Box of 10

Stock No
7331

Hollister Universal Catheter Access Port



9779

- ▶ Graduated nipple and punch allows introduction of catheter or drain through pouch
- ▶ Sterile

Box of 8

Stock No
9779

Other Accessories

Loop Ostomy Bridge



7767

- ▶ Sterile

HCPCS: A4421

Box of 20

Stock No

7767

Bedside Drainage Collection System with Anti-Reflux Valve



9839

- ▶ Sample Port
- ▶ Sterile Fluid Pathway

HCPCS: A4357

Box of 20

Stock No

Description

9839

68 oz (2000 mL) Capacity

Skin Barrier Pastes and Powders



7910

- ▶ Karaya 5 Paste

HCPCS: A4405

Box of 12

Stock No

Description

7910

4.5 oz (135 mL) tube



7905

- ▶ Karaya 5 Powder

HCPCS: A4371

Box of 12

Stock No

Description

7905

2.5 oz (75 mL) puff bottle

Skin Barriers



7800

- ▶ Premium (standard wear) Skin Barrier

HCPCS: A4362

Box of 5

Stock No

Size

7800

4 x 4" (10 cm x 10 cm)



7700

- ▶ HolliHesive (standard wear) Skin Barrier

HCPCS: A4362

Box of 5

Stock No

Size

7700

4 x 4" (10 cm x 10 cm)

Replacement Filters



7766

- ▶ For use with 3668
- ▶ 1 pack contains 20 Filters and 2 Plugs

HCPCS: A4368

Box of 5 Packs

Stock No

7766

Vented Urinary Leg Bag Combination Pack



9645/9655

Box of 10

► Includes:

- (1) Sterile, Vented Urinary Leg Bag — Medium or Large
- (1) Sterile, Extension Tubing with Vacuum Relief Valve, (18" Long, Oval, 1 1/32" Effective Inside Diameter), with Connector
- (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)

► Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length	Width
9645	Medium — 18 oz (540 mL)	10"	4-1/4"
9655	Large — 30 oz (900 mL)	11-1/8"	5"

Urinary Leg Bag Combination Pack



9348/9349

Box of 12

► Includes:

- (1) Non-Sterile, Urinary leg Bag — Medium or Large
- (1) Non-Sterile, Extension Tubing (18" Long, Oval, 11/32" Effective Inside Diameter) with Connector
- (1) Pair, Fabric Leg Bag Straps — Medium (23" Long, 1-1/8" Wide)

► Not Made with Natural Rubber Latex

HCPCS: A4358

Stock No	Description	Length	Width
9348	Medium — 18 oz (540 mL)	10"	4-1/4"
9349	Large — 30 oz (900 mL)	11-1/8"	5"

Urinary Leg Bag



9814/9805

Box of 10

► Includes:

- (1) Sterile, Urinary Leg Bag — Medium or Large

► Not Made with Natural Rubber Latex

HCPCS: A9270

HCPCS: A4358 (if straps are added)

Stock No	Description	Length	Width
9814	Medium — 18 oz (540 mL)	10"	4-1/4"
9805	Large — 30 oz (900 mL)	11-1/8"	5"

Fabric Leg Bag Straps



9342/9343

10 Pair

- ▶ Not Made with Natural Rubber Latex

HCPCS: A5114

Stock No	Description	Length	Width
9342	Small	15"	1-1/8"
9343	Medium	23"	1-1/8"

Extension Tubing



9345/9346

Box of 10


- ▶ Kink-resistant with Connector
- ▶ Oval, 11/32" Effective Inside Diameter
- ▶ Thermoplastic Elastomer
- ▶ Not Made with Natural Rubber Latex

HCPCS: A4331

Stock No	Description	Length
9345	Non Sterile	18"
9346	Sterile	18"

General Information

- ▶ Similar to an enema with the water being placed into the stoma.
- ▶ Routine irrigation is recommended to be done on doctor's orders (to regulate emptying the bowel at a scheduled time).
- ▶ Appropriate for some people with sigmoid or descending colostomies.
- ▶ Not appropriate for people with urostomy or ileostomy.
- ▶ Allow approximately one hour for irrigation procedure.



7721

- ▶ Cone Irrigator Kit Includes:
 - ▶ (1) Irrigator
 - ▶ (1) Stoma Cone with Connector
 - ▶ (1) Cleaning Brush

Box of 1

**Irrigation Bag Holds
Approximately
1500 CC**

HCPCS: A4398 + A4399

Stock No

7721



7723

- ▶ Stoma Cone with Connector
- ▶ Replacement Unit

Box of 10

HCPCS: A4399

Stock No

7723



1460x (CUT-TO-FIT)

- ▶ *Flexend* (extended wear) Skin Barrier
- ▶ With Tape Border




Box of 5

HCPCS: A4409

Stock No

Flange Size

For Stomas Up To

14602	 1-3/4" (44 mm)	1-1/4" (32 mm)
14603	 2-1/4" (57 mm)	1-3/4" (44 mm)
14604	 2-3/4" (70 mm)	2-1/4" (57 mm)



1815x

- ▶ Drainable Irrigator Sleeve to be used with *New Image* Two-Piece System
- ▶ Fits any New Image barrier that has a 1-3/4", 2-1/4" or 2-3/4" flange, including the 1460x series above
- ▶ With Belt Tabs
- ▶ Includes Clamp

Box of 5

HCPCS: A4397

Transparent

Stock No

Flange Size

Sleeve Length

18152	 1-3/4" (44 mm)	30" (76 cm)
18153	 2-1/4" (57 mm)	30" (76 cm)
18154	 2-3/4" (70 mm)	30" (76 cm)

**7728/7724****Box of 20**

- ▶ Irrigator Sleeve with Belt Tabs
- ▶ Includes Clamp

HCPCS: A4397

Stock No	Opening Size	Sleeve Length
7728	2" (51 mm)	22" (56 cm)
7724	3" (76 mm)	22" (56 cm)

**318x (PRE-SIZED)****Box of 30**

- ▶ Stoma Cap, Tape Adhesive, Integrated Filter, Odor-Barrier Film, Without *ComfortWear* Panels

HCPCS: A5055

Beige Stock No	Opening	Size
1796	1-15/16" (50 mm)	4" (102 mm)
3184	2" (51 mm)	4-1/4" (11 cm)
3186	3" (76 mm)	4-1/4" (11 cm)

**7300/7299****Box of 10**

- ▶ Ostomy Belts

HCPCS: A4367

Beige Stock No	Size
7300	23"-43" (58 cm-109 cm) Medium
7299	34"-65" (86 cm-165 cm) Large

**7740****Box of 10**

- ▶ *Adapt* Stoma Lubricant
- ▶ For Use with Stoma Irrigator

HCPCS: A4402

Stock No	Description
7740	4 oz (120 mL) bottle

For a complete list of Hollister products, refer to the comprehensive Hollister product catalog online at www.hollister.com.

Karaya 5 One-Piece Pouching System

Drainable Pouches



366x (PRE-SIZED)

- ▶ Drainable Pouch, *Karaya 5* (standard wear) Skin Barrier — Tape Border, Replaceable Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels. Includes 1 Bag of 20 Hollister Replacement Filters and 2 Plugs, #7766

Box of 10



HCPCS: A4389 + A4368

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3668	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)



360x (PRE-SIZED)

- ▶ Drainable Pouch, *Karaya 5* (standard wear) Skin Barrier — Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

Box of 10



HCPCS: A4389

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3608	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)
3603	1-1/2" (38 mm)	1-1/8" (29 mm)	12" (30 cm)
3609	1-3/4" (44 mm)	1-3/8" (35 mm)	12" (30 cm)
3604	2" (51 mm)	1-1/2" (38 mm)	12" (30 cm)
3605	2-1/2" (64 mm)	2" (51 mm)	12" (30 cm)



323x (PRE-SIZED)

- ▶ Drainable Mini-Pouch, *Karaya 5* (standard wear) Skin Barrier — Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

Box of 30



HCPCS: A4389

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3232	1" (25 mm)	5/8" (16 mm)	9" (23 cm)
3238	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
3233	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
3239	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
3234	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)



322x (PRE-SIZED)

- ▶ Drainable Pouch, *Karaya 5* (standard wear) Skin Barrier — Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

Box of 30



HCPCS: A4389

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3222	1" (25 mm)	5/8" (16 mm)	12" (30 cm)
3228	1-1/4" (32 mm)	7/8" (22 mm)	12" (30 cm)
3223	1-1/2" (38 mm)	1-1/8" (29 mm)	12" (30 cm)
3229	1-3/4" (44 mm)	1-3/8" (35 mm)	12" (30 cm)
3224	2" (51 mm)	1-1/2" (38 mm)	12" (30 cm)
3225	2-1/2" (64 mm)	2" (51 mm)	12" (30 cm)
3226	3" (76 mm)	2-1/2" (64 mm)	12" (30 cm)

Drainable Pouches (continued)



327x (PRE-SIZED)

- ▶ Drainable Pouch, *Karaya 5* (standard wear)
Skin Barrier — Tape Border, Without Integrated Filter, Pouch Clamp Closure, Belt Tabs, Without *ComfortWear* Panels

Box of 30



HCPCS: A4389

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3278	1-1/4" (32 mm)	7/8" (22 mm)	16" (40 cm)
3273	1-1/2" (38 mm)	1-1/8" (29 mm)	16" (40 cm)
3279	1-3/4" (44 mm)	1-3/8" (35 mm)	16" (40 cm)
3274	2" (51 mm)	1-1/2" (38 mm)	16" (40 cm)
3275	2-1/2" (64 mm)	2" (51 mm)	16" (40 cm)
3276	3" (76 mm)	2-1/2" (64 mm)	16" (40 cm)

Closed Pouches



332x (PRE-SIZED)

- ▶ Closed Pouch, *Karaya 5* (standard wear)
Skin Barrier — Tape Border, Integrated Filter, Belt Tabs, Without *ComfortWear* Panels

Box of 30



HCPCS: A4417

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
3328	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
3323	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
3329	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
3324	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)
3325	2-1/2" (64 mm)	2" (51 mm)	9" (23 cm)
3326	3" (76 mm)	2-1/2" (64 mm)	9" (23 cm)



716x (PRE-SIZED)

- ▶ Closed Pouch, *Karaya 5* (standard wear)
Skin Barrier — Without Tape Border, Without Integrated Filter, Belt Tabs, Without *ComfortWear* Panels

Box of 30



- ▶ Requires a Belt

HCPCS: A4387

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
7164	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)

Urostomy Pouches



143x (PRE-SIZED)

- ▶ Urostomy Pouch, *Karaya 5* (standard wear)
Skin Barrier — Tape Border, Belt Tabs, Without *ComfortWear* Panels

Box of 10



HCPCS: A4429

Transparent Stock No	Pouch Size	Barrier Opening	Pouch Length
1432	1" (25 mm)	5/8" (16 mm)	9" (23 cm)
1438	1-1/4" (32 mm)	7/8" (22 mm)	9" (23 cm)
1433	1-1/2" (38 mm)	1-1/8" (29 mm)	9" (23 cm)
1439	1-3/4" (44 mm)	1-3/8" (35 mm)	9" (23 cm)
1434	2" (51 mm)	1-1/2" (38 mm)	9" (23 cm)

A young woman with long brown hair, wearing a denim jacket, is sitting at a wooden table and looking at a silver tablet. An older man with a grey beard and a grey sweater is sitting next to her, also looking at the tablet. They appear to be in a bright, indoor setting with large windows in the background showing greenery.

Peristomal skin should be

as healthy as the skin on

the rest of the body

Education

Education

All of the Hollister educational materials listed below are also available on our website: www.hollister.com



Visit the Hollister YouTube channel for all educational videos

Case Studies

For a complete description and current listing, visit the Clinical Evaluations and Case Studies section found in the Professional Resources portion of our website, www.hollister.com. For more information, contact your Hollister Ostomy Sales Representative.

Stoma Measuring Guide	903515
Oval Stoma Measuring Guide	923298
<i>(Printable versions available online at www.hollister.com/en/ostomycare/educationaltools)</i>	

Educational Booklets	Literature No
Understanding Your Colostomy	923054
Understanding Your Colostomy (Spanish) <i>(online only)</i>	923234
Understanding Your Ileostomy	923055
Understanding Your Ileostomy (Spanish) <i>(online only)</i>	923235
Understanding Your Urostomy	923056
Understanding Your Urostomy (Spanish) <i>(online only)</i>	923236
Can They Still Wiggle and Giggle? Coloring Book	907617
The Day I Got My Stoma <i>(online only)</i>	910497
Caring for Your Child With an Ostomy – Pediatric (English)	923057
Caring for Your Child With an Ostomy – Pediatric (Spanish) <i>(online only)</i>	923257
Caring for your Loved One with an Ostomy	923058
Caring for your Loved One with an Ostomy (Spanish) <i>(online only)</i>	923258




DID YOU KNOW...

- All of the educational resources listed are also available online at www.hollister.com.
- Your patients can also access a library of educational articles on topics related to ostomy care - ranging from general ostomy information and maintaining healthy peristomal skin, to lifestyle changes and how to adapt to life with a stoma.



Ostomy Care Tips	Literature No
<i>Adapt</i> Barrier Rings	923059
<i>Adapt</i> Convex Barrier Rings	923060
<i>Adapt</i> Lubricating Deodorant	923061
<i>Adapt</i> Medical Adhesive Spray/ <i>Adapt</i> Remover Spray & Wipes	923073
<i>Adapt</i> No-Sting Skin Protective Wipe	923079
<i>Adapt</i> Ostomy Belt	923062
<i>Adapt</i> Paste	923063
<i>Adapt</i> Stoma Powder	923064
Colostomy Irrigation (<i>online only</i>)	923065
Using Convexity	923066
Fistula Management (<i>online only</i>)	923067
High Output	923121
Ileostomy Dietary Guidelines	923117
<i>Lock 'n Roll</i> Microseal Closure	923068
Medicare Reimbursement Guide	923118
One-Piece Pouching System (English)	923069
One-Piece Pouching System (Spanish)	923070
Ostomy Educational Theatre (<i>online only</i>)	923115
Patient Anatomy Guide	923072
Peristomal Skin Care	923085
<i>Pouchkins</i> Newborn Pouching System (<i>online only</i>)	923119
<i>Pouchkins</i> One-Piece Pouching System (<i>online only</i>)	923074
<i>Pouchkins</i> Premie Pouching System (<i>online only</i>)	923076
<i>Pouchkins</i> Two-Piece Pouching System (<i>online only</i>)	923075
Routine Care of Your Ostomy (English)	923077
Routine Care of Your Ostomy (Spanish)	923078
Stoma Site Selection (<i>online only</i>)	923080
Two-Piece Pouching System (English)	923081
Two-Piece Pouching System (Spanish)	923082
Two-Piece Pouching System with <i>FormaFlex</i> Skin Barrier (<i>online only</i>)	923083
Types of Ostomies	923084
Universal Catheter Access Port (UCAP) (<i>online only</i>)	923114
Where a Soft Convex Skin Barrier May Improve Fit	923116
Spanish Lit QR sheet	923398
Patient QR Code Sheet	923423
Educational Theater Videos	923429

Living with an Ostomy Booklets	Literature No
Travel	923123
Healthy Eating	923124
Sex & Parenthood	923125
Home & Work Life	923126
Sports & Fitness	923127
<i>Adapt</i> Barrier Extender	923899

Pre-Surgery Care Kits	Literature No
	<p>Includes guidebook to facilitate the patient's pre-surgery visit with a clinician and materials for self-led practice applying and removing an ostomy pouching system with a model stoma.</p>
Colostomy Pre-Surgery Care Kit	K2025
Ileostomy Pre-Surgery Care Kit	K2026
Urostomy Pre-Surgery Care Kit	K2027

Patient Discharge Folders

Colostomy Patient Discharge Folder	923261
Ileostomy Patient Discharge Folder	923262
Urostomy Patient Discharge Folder	923263

Professional Organizations

Wound, Ostomy and Continence Nurses Society (WOCN)
www.wocn.org
1.888.224.9626

Professional Resources

Fit Indicator Tool (FIT)

Fit Indicator Tool (FIT) has been designed to help provide a simple patient assessment methodology to determine the need for a flat or convex barrier. Open the app Store or Google Play and search using: "Hollister Education".

Peristomal Skin Assessment Guide

The Peristomal Skin Assessment Guide helps with assessment and management of common peristomal skin complications and is available at the WOCN Society website. Also available for consumer use. psag.wocn.org

Online Education Program

Our Ostomy Care course allows healthcare professionals to learn about different types of ostomies, problem-solving tips, and additional resources. It also provides flexibility and the ease of online, anytime e-learning.
www.hollister.com/ostomyelearning

Shadow Buddies



	Item No
Light Skin Tone Girl	905439
Light Skin Tone Boy	905440
Medium Skin Tone Girl	905441
Medium Skin Tone Boy	905442
Dark Skin Tone Girl	905443
Dark Skin Tone Boy	905444

Additional Resources

	Item No
Stoma Model (White Plastic)	858674

Consumer Resources

Hollister Secure Start Services
1.888.808.7456
www.securestartservices.com

United Ostomy Associations of America, Inc. (UOAA)
info@uoaa.org
www.ostomy.org
1.800.826.0826

Pull-Thru Network
pullthrunetwork.org
1-309-262-0786

Crohn's and Colitis Foundation of America (CCFA)
www.cdfa.org
1.800.932.2423

Shadow Buddies
www.shadowbuddies.org
1.913.642.4646

AF300 Filter: “AF” stands for air flow and “300” represents the 300 mm² air flow surface area on every AF300 filter.

Ceramide: A natural component of the human skin that links the cells of the stratum corneum together to form a waterproof protective barrier. This helps prevent water loss that can lead to skin damage and dryness.

CeraPlus Skin Barrier: A skin barrier that has all the basic characteristics of an ostomy skin barrier (adhesion, absorption, and erosion resistance) and is infused with ceramide. *CeraPlus* is an example of an extended wear skin barrier.

Closed Pouch: A pouch without an opening at the bottom. Quick and easy to discard. It must be removed to be emptied.

Colostomy: An ostomy (surgical opening) created in the colon, part of the large intestine.

ComfortWear (pouch) Panels: Provide a soft cloth between the pouch film and the skin.

Convex Skin Barrier: A skin barrier shaped to press more deeply around the stoma. May be used when the stoma does not protrude from the abdomen.

Cut-to-Fit Skin Barrier: A skin barrier that can be customized to a needed size and shape. Ideal for stomas that are not round or stomas that are still changing in size.

Drainable Pouch: A pouch with an opening at the bottom. An integrated closure system or plastic clamp is used to keep the pouch closed until it is time to empty it.

Filter: Releases gas but not odor out of the pouch. The key performance factors of any filter are air flow, deodorization, and liquid protection.

Flange: The plastic ring in a two-piece ostomy system that is used to connect the skin barrier and pouch together.

Flat Skin Barrier: May be used when the stoma protrudes from the abdomen.

Flexend Skin Barrier: A skin barrier designed to be durable and resistant to liquid discharge. *Flexend* is an example of an extended wear skin barrier.

FlexWear Skin Barrier: A skin barrier designed to be worn more than one day. A *FlexWear* skin barrier is designed to be more resistant to breakdown than a *SoftFlex* skin barrier. *FlexWear* is an example of a standard wear skin barrier.

Floating Flange: Type of skin barrier flange found on Hollister two-piece pouching systems that allows you to put your fingers under the skin barrier flange while attaching the pouch.

FormaFlex Skin Barrier: A skin barrier that can be stretched to fit the size and shape of the stoma without the use of scissors. It is a barrier option in the New Image two-piece pouching system. *FormaFlex* is an example of an extended wear skin barrier.

Ileostomy: A surgically created opening in the last portion of the small intestine.

Irrigation: A procedure sometimes used by individuals with a colostomy to control bowel elimination. It is not for everyone with a colostomy and requires special equipment and instruction by a knowledgeable healthcare provider.

Lock 'n Roll Microseal Closure: Integrated system for closing the bottom of drainable pouches. Utilizes plastic interlocking fasteners to create a secure closure. Found on *New Image*, *Premier*, and *Pouchkins* drainable pouches.

Medical Adhesive-Related Skin Injury (MARS): Occurrences in which erythema and/or other manifestation such as blister, tear or other injury persists more than 30 minutes after removal of medical adhesive.

Odor-Barrier Pouch: Is highly resistant to letting odors pass through the pouch film. The odor barrier quality increases confidence when wearing a pouch.

Peristomal Moisture-Associated Skin Damage (PMASD): Inflammation and erosion of the skin adjacent to the stoma, associated with exposure to effluent such as urine or stool.

Pre-Sized Skin Barrier: A skin barrier that is already cut to the proper size.

QuietWear Pouch Material: Material comprised of a pouch film and *ComfortWear* panel that are laminated together. This special design results in the reduction of the pouch rustling noise.

SoftFlex Skin Barrier: A skin barrier that provides gentle adhesive attachment to the skin. *SoftFlex* is designed to limit skin damage with frequent removal and can be changed daily (or more often). *SoftFlex* is an example of a standard wear skin barrier.

Stoma: A surgically created opening in the gastrointestinal or urinary tract. Also known as an ostomy. The stoma is pink and moist and doesn't have feeling.

Urostomy: A surgically created opening into the urinary system; drains urine and sometimes mucus.

Urostomy Pouches: Pouch with a drain tap at the bottom so urine can be emptied quickly and easily.

Wear Time: The length of time a product can be worn before failure. Wear time varies widely based upon a variety of factors.

Key Goals for Nursing Management

- **Protect the skin** – The output from a fistula can be very irritating to the skin. Broken skin contributes to pain and discomfort. A skin barrier can protect the skin from contact with the drainage.
- **Contain the output** – If the output is minimal, gauze dressings may be sufficient for containment. However, if the volume of the output from the fistula is significant or frequently overwhelms the dressing, a pouch may be more appropriate. Using a pouch can be more efficient and can protect healthcare providers from exposure to bodily fluids.
- **Measure the output** – With a high output fistula, loss of fluids and electrolytes can be a major concern. Measuring the drainage is easier using a pouching system. Sometimes connecting a pouch to a bedside collector can relieve the stress on the pouch, so the wear time improves.
- **Manage costs** – Care of fistulas can be costly in terms of supplies and human resources. Using a pouching system and appropriate accessories can be more cost-effective and manageable than frequent dressing changes.
- **Control odor** – Containment of the drainage in a pouching system, even in small quantities, can manage odor effectively. Other accessory products can also assist with odor control when the pouch is changed or emptied.
- **Provide patient comfort** – Effectively managing a fistula can contribute to patient dignity. Protecting the skin, collecting the discharge, containing costs, and controlling odor all contribute to patient comfort.

Index • Stock Number

<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>
1120x	15, 29	18282	38, 75	366x	88
1140x	15, 32	1829x	38	3761	75
1150x	15, 34	18292	38, 75	3774	76
1340x	33	1832x	40	3777	73
1350x	35	1833x	41	3778	73
139xx	34	1835x	41	3795	74
1410x	31	1836x	40	3796	73
1420x	29	1837x	40	3797	74
143x	89	1838x	41	3799	75
1430x	31	18382	41, 75	7164	89
1440x	33	1839x	41	7299	80
1450x	35	18392	41, 75	7300	80
1460x	86	1840x	42	7331	82
1470x	30	18402	42, 76	7700	83
1480x	32	1841x	42	771x	82
149xx	34	18412	42, 76	7721	86
1510x	15, 29	1842x	42	7723	86
1520x	30	18422	42, 76	7724	87
1530x	32	1873x	41	7728	87
1540x	33	1875x	41	773x	82
1550x	34	1890x	42	7730	81
1560x	29	1891x	42	7736	82
1580x	32	1892x	42	7737	81
1590x	35	1900x	43	7740	87
1610x	30	1905x	45	7760	81
1640x	31	1910x	43	7766	83
1796	65, 87	1915x	45	7767	83
1800x	38	1920x	44	78xx	78
1801x	40	1925x	46	7800	83
1810x	39	1930x	44	78500	80
1811x	38	1935x	45	78501	80
1812x	39	1940x	44	7905	83
1813x	38	1950x	46	7906	80
1814x	39	1960x	21, 43	7910	83
1815x	86	1980x	21, 43	7917	81
1816x	39	318x	87	79300	80
1817x	39	322x	88	79301	80
1818x	38	323x	88	79400	80
1819x	38	327x	89	79402	80
1820x	39	332x	89	795xx	78
1828x	38	360x	88	796xx	78

<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>	<u>Stock No</u>	<u>Page(s)</u>
80070	63	849x	69	89003	71
80110	63	8499x	69	89004	71
8031	50	850xx	58	89005	71
8081	49	851xx	61	891x	19, 56
8131	50	852xx	58	891xx	19, 56
8181	49	853x	52	8931	17, 48
821xx	65	8531	48	895xx	18, 23, 55, 78
82100	64	8/54x	51	896x	18, 55
823xx	64	855x	52	896xx	23, 78
823x00	58	8551	48	898x	19, 56
82302	64	857xx	55	898xx	19, 56
824xx	58	858xx	57	899x	19, 60
82400	64	859x	61	970x	63
82402	64	861x	62	977x	63
82500	64	862xx	59	9779	82
8261x	62	8631	49	9839	83
8262x	59	863x	52		
827xx	57	8641	49		
828x	61	864x	52		
829x	61	8651	49		
8331	48	866x	56		
8371	48	867xx	55		
8381	49	8770	82		
841x	20, 70	88xx	23, 78		
8413xx	20, 70	8800	81		
843xx	68	8801	81		
8439xx	68	883xx	53		
844x	67	88300	50		
8440	66	88302	51		
8449x	67	884xx	53		
84490	66	88400	50		
8450	66	88401	51		
84590	66	88402	51		
8460	66	88501	51		
846x	67	887xx	53		
8469x	67	88700	50		
84690	66	88800	50		
847xx	68	89xx	17, 51		
8479xx	68	8901	17, 48		
848x	69	89001	71		
8489x	69	89002	71		

Name	Page(s)
A dapt Barrier Extenders	80
<i>Adapt</i> CeraRing Convex Barrier Rings.....	23, 78
<i>Adapt</i> CeraRing Flat Barrier Rings	23, 78
<i>Adapt</i> CeraRing Oval Convex Barrier Rings	23, 78
<i>Adapt</i> Odor Eliminator and Lubricant.....	80
<i>Adapt</i> Ostomy Belts	80
<i>Adapt</i> Skin Barrier Rings.....	78
<i>Adapt</i> Skin Barrier Paste.....	80
<i>Adapt</i> Skin Barrier Strips.....	80
<i>Adapt</i> Skin Prepping Product	81
<i>Adapt</i> Stoma Powder.....	80
<i>Adapt</i> Adhesives and Removers.....	81
B edside Drainage Collection Systwm with Anti-Reflux Valve	83
C eraPlus Convex Skin Barriers with Integrated Floating Flange	15-16, 32, 34
<i>CeraPlus</i> Soft Convex Skin Barriers with Integrated Floating Flange.....	16, 36
<i>CeraPlus</i> Flat Skin Barriers with Integrated Floating Flange	15, 29
E xtension Tubing	85
D rainable Pouch Clamp	82
Drainable Pouches with <i>CeraPlus</i> Flat Barriers.....	17, 48, 51
Drainable Pouches with <i>CeraPlus</i> Soft Convex Barriers	18, 55-56
Drainable Pouches with <i>CeraPlus</i> Convex Barriers.....	19, 56, 60
F lex <i>extend</i> (extended wear) Skin Barrier.....	81
Fabric Leg Bag Straps	85
H igh Output Drainable Pouches.....	40, 63
<i>Hollister</i> Universal Catheter Access Port	82
I rrigation Products.....	86-87
Irrigator Sleeves	86-87
K araya 5 Closed Pouches	89
<i>Karaya</i> 5 Drainabl Pouches.....	88-89
<i>Karaya</i> 5 Urosotmy Pouches	89
L oop Ostomy Bridge	83
m 9 Cleaner/Decrystallizer	82
<i>m9</i> Odor Eliminators.....	82

N <i>New Image</i> Closed Pouches.....	40-41
<i>New Image</i> Soft Convex Skin Barrier with Integrated Floating Flange.....	36
<i>New Image</i> Convex Skin Barriers with Integrated Floating Flange.....	32-35
<i>New Image</i> Drainable Pouches.....	38-39
<i>New Image</i> Flat Skin Barriers with Integrated Floating Flange.....	29-31
<i>New Image</i> High Output Drainable Pouches.....	40
<i>New Image</i> Non Sterile Kits - Drainable.....	43-44
<i>New Image</i> Non Sterile Kits - Urostomy.....	44
<i>New Image</i> Sterile Kits - Drainable.....	45
<i>New Image</i> Sterile Kits - Urosotmy.....	46
<i>New Image</i> Urostomy Pouches.....	42
P <i>ouchkins</i> Newborn Pouch.....	73
<i>Pouchkins</i> One-Piece Pouches.....	73-74
<i>Pouchkins</i> Pediatric Ostomy Belt.....	76
<i>Pouchkins</i> Premie Pouch.....	73
<i>Pouchkins</i> Two-Piece Pouches.....	75-76
<i>Pouchkins</i> Two-Piece Skin Barrier.....	75
<i>Premier</i> Closed Pouches with Flat Barriers.....	64-65
<i>Premier</i> Drainable Pouches with Convex Barriers.....	55-62
<i>Premier</i> Drainable Pouches with Flat Barriers.....	48-53
<i>Premier</i> High Output Pouches with Soft Convex Barriers.....	55-56
<i>Premier</i> High Output Drianable Pouches.....	63
<i>Premier</i> Non Sterile Kits - Drainable.....	71
<i>Premier</i> Non Sterile Kits - Urosotmy.....	71
<i>Premier</i> Urostomy Pouches with Convex Barriers.....	68-70
<i>Premier</i> Urosotmy Pouches with Flat Barriers.....	66-67
R eplacement Filters.....	83
S kin Barriers.....	83
Skin Barrier Pastes and Powders.....	83
Stoma Cap.....	65, 87
U rostomy Drain Tube Adapter.....	82
Urosotmy Pouches with <i>CeraPlus</i> Soft Convex Barriers.....	20, 70
Urinary Leg Bags.....	84
W ound Drainage Collectors.....	63

General Information

HCPCS Codes Disclaimer

The reimbursement information provided herein is intended to provide general information concerning coding of Hollister products only. Hollister does not guarantee coverage or payment for any products. The ultimate responsibility for proper coding, satisfying reimbursement requirements, and obtaining reimbursement remains with the provider. Coding and coverage policies and guidelines are complex, can vary from one carrier or region to another, and are updated frequently. Providers should check with their local carriers or intermediaries often and should consult with counsel, a reimbursement specialist, and/or the DMECS website for any coding, coverage, reimbursement, or billing questions.

HCPCS codes shown in this Product Catalog are current as of January 2018.

General Medicare Product Maximum Allowables of Selected Ostomy Products

Two-Piece Pouching Systems	Maximum Quantity per Month
Flat, Standard Wear Skin Barriers with Floating Flanges	20
Drainable Pouches	20
Closed Pouches	60
Urostomy Pouches	20
Irrigator Sleeves	4

One-Piece Pouching Systems	Maximum Quantity per Month
Drainable Pouches with Flat, Standard Wear Skin Barriers	20
Closed Pouches with Flat, Standard Wear Skin Barriers	60
Stoma Caps	31

Pediatric Pouching Systems	Maximum Quantity per Month
One-Piece Drainable and Urostomy Pouches	20
Two-Piece Skin Barriers	20
Two-Piece Drainable and Urostomy Pouches	20
Two-Piece Mini Closed Pouches	60
Pediatric Ostomy Belt	1

Accessories	Maximum Quantity per Month
Skin Barrier Rings	*
Skin Barrier Paste	4 oz
Stoma Powder	10 oz/6 months
Lubricating Deodorant	*
Ostomy Belt	1
Skin Protective Wipes	150/6 months

Note: Medicare allows a usual maximum quantity of supplies each month. The table above provides the maximum allowable amounts. The amount of supplies you need may vary from the maximum quantities. That will be determined by your physician and must be on your order and in your medical record. If you require more ostomy supplies than what is listed in the table, your physician will need to provide a medical justification for the additional supplies. The ostomy supplier will need to have that medical justification on file when billing the Medicare program for the additional supplies.

* No usual maximum amount established.

Reimbursement For Fistula Management by Medicare

Reimbursement for fistula management under the Medicare program is described in the Surgical Dressings Policy Article. Specifically the article states:

“Surgical dressings are covered when either of the following criteria are met:

1. They are required for the treatment of a wound caused by, or treated by, a surgical procedure; **or**
2. They are required after debridement of a wound.”

The article goes on to explain:

“Examples of situations in which dressings are noncovered under the Surgical Dressings benefit are:

- a. Drainage from a cutaneous fistula which has not been caused by or treated by a surgical procedure”

For information see The Fistula Management Care Tip on www.hollister.com

Latex Statement

All Hollister Ostomy Care products found in this catalog and its packaging does not include natural rubber latex as a component, and our component suppliers do not add natural rubber as part of their production process.

Prior to use, be sure to read the Instructions for Use for information regarding intended use, contraindications, warnings, precautions and instructions.



DID YOU KNOW...

Supplier Matching

“Having a supplier that can provide your preferred products when you need them is important. If you are having a difficult time finding a supplier who provides your preferred products, please call 888-808-7456 to learn your options.”

– Secure Start Member

Hollister Secure Start services can connect patients with a supplier of their choice based on their insurance plan for the products you recommend. The amount of supplies they need may vary from the maximum allowable quantities, but this is determined by a physician and should be included on the patient's order and medical record. For additional assistance, visit www.securestartservices.com or call us at 1.888.808.7456.



Medicare's Usual Maximum Quantity for Ostomy Supplies

The quantity of ostomy supplies needed by a patient is determined primarily by the type of ostomy, its location, its construction, and the condition of the skin surface surrounding the stoma. There will be variation according to individual patient need and their needs may vary over time. The table below lists the maximum number of items/units of service that are usually medically necessary. The actual quantity needed for a particular patient may be more or less than the amount listed depending on the factors that affect the frequency of barrier and pouch change. The medical necessity for use of a greater quantity of supplies than the amounts listed must be clearly documented in the patient's medical record, and must be available upon request. If adequate documentation is not provided when requested, the excess quantities will be denied as not medically necessary.

Per Medicare LCD for Ostomy Supplies (L33828), effective 1-1-2020.

HCPCS	Long Description	Maximum Quantity
A4331	Extension drainage tubing, any type, any length, with connector/adapter, for use with urinary leg bag or urostomy pouch, each	*
A4357	Bedside drainage bag, day or night, with or without anti-reflux device, with or without tube, each	2 ea /mo
A4361	Ostomy faceplate, each	3 ea/6 mo
A4362	Skin barrier, solid, 4 x 4 or equivalent, each	20 ea/mo
A4363	Ostomy clamp, any type, replacement only, each	*
A4364	Adhesive, liquid or equal, any type, per ounce	4 oz/mo
A4366	Ostomy vent, any type, each	*
A4367	Ostomy belt, each	1 ea/mo
A4368	Ostomy filter, any type, each	*
A4369	Ostomy skin barrier, liquid (spray, brush, etc.), per ounce	2 oz/mo
A4371	Ostomy skin barrier, powder, per ounce	10 oz/6 mo
A4372	Ostomy skin barrier, solid 4 x 4 or equivalent, standard wear with built-in convexity, each	*
A4373	Ostomy skin barrier, with flange (solid, flexible, or accordion), with built-in convexity, any size, each	*
A4375	Ostomy pouch, drainable, with faceplate attached, plastic, each	*
A4376	Ostomy pouch, drainable, with faceplate attached, rubber, each	*
A4377	Ostomy pouch drainable for use on faceplate, plastic, each	10 ea/mo
A4378	Ostomy pouch, drainable, for use on faceplate, rubber, each	*
A4379	Ostomy pouch, urinary, with faceplate attached, plastic, each	*
A4380	Ostomy pouch, urinary, with faceplate attached, rubber, each	*
A4381	Ostomy pouch, urinary, for use with faceplate, plastic, each	10 ea/mo
A4382	Ostomy pouch, urinary, for use on faceplate, heavy plastic, each	*
A4383	Ostomy pouch, urinary, for use on faceplate, rubber, each	*
A4384	Ostomy faceplate equivalent, silicone ring, each	*
A4385	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, without built-in convexity, each	*
A4387	Ostomy pouch, closed, with barrier attached, with built-in convexity (one-piece), each	*
A4388	Ostomy pouch, drainable, with extended wear barrier attached, (one-piece), each	*
A4389	Ostomy pouch, drainable, with barrier attached, with built-in convexity (one-piece), each	*

*No established usual maximum, refer to Medicare Utilization Guidelines.

HCPGS	Long Description	Quantity
A4390	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity (one-piece), each	*
A4391	Ostomy pouch, urinary, with extended wear barrier attached (one-piece), each	*
A4392	Ostomy pouch, urinary, with standard wear barrier attached, with built-in convexity (one-piece), each	*
A4393	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity (one-piece), each	*
A4394	Ostomy deodorant, with or without lubricant, for use in ostomy pouch, liquid, per fluid ounce	*
A4395	Ostomy deodorant for use in ostomy pouch, solid, per tablet	*
A4396	Ostomy belt with peristomal hernia support	*
A4397	Irrigation supply, sleeve, each	4 ea/mo
A4398	Ostomy irrigation supply, bag, each	2 ea/6 mo
A4399	Ostomy irrigation supply, cone/catheter, with or without brush	2 ea/6 mo
A4402	Lubricant, per ounce	4 oz/mo
A4404	Ostomy ring, each	10 ea/mo
A4405	Ostomy skin barrier, non-pectin based, paste, per ounce	4 oz/mo
A4406	Ostomy skin barrier, pectin-based, paste, per ounce	4 oz/mo
A4407	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, 4 x 4 inches or smaller, each	*
A4408	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, with built-in convexity, larger than 4 x 4 inches, each	*
A4409	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, without built-in convexity, 4 x 4 inches or smaller, each	*
A4410	Ostomy skin barrier, with flange (solid, flexible, or accordion), extended wear, without built-in convexity, larger than 4 x 4 inches, each	*
A4411	Ostomy skin barrier, solid 4 x 4 or equivalent, extended wear, with built-in convexity, each	*
A4412	Ostomy pouch, drainable, high output, for use on a barrier with flange (two-piece system), without filter, each	*
A4413	Ostomy pouch, drainable, high output, for use on a barrier with flange (two-piece system), with filter, each	*
A4414	Ostomy skin barrier, with flange (solid, flexible, or accordion), without built-in convexity, 4 x 4 inches or smaller, each	20 ea/mo
A4415	Ostomy skin barrier, with flange (solid, flexible, or accordion), without built-in convexity, larger than 4 x 4 inches, each	20 ea/mo
A4416	Ostomy pouch, closed, with barrier attached, with filter (one-piece), each	60 ea/mo
A4417	Ostomy pouch, closed, with barrier attached, with built-in convexity, with filter (one-piece), each	60 ea/mo
A4418	Ostomy pouch, closed, without barrier attached, with filter (one-piece), each	60 ea/mo
A4419	Ostomy pouch, closed, for use on barrier with non-locking flange, with filter (two-piece), each	60 ea/mo
A4420	Ostomy pouch, closed, for use on barrier with locking flange (two-piece), each	60 ea/mo
A4421	Ostomy supply; miscellaneous	*
A4422	Ostomy absorbent material (sheet/pad/crystal packet) for use in ostomy pouch to thicken liquid stomal output, each	*

* No established usual maximum, refer to Medicare Utilization Guidelines.

Medicare Utilization Chart

HCPCS	Long Description	Quantity
A4423	Ostomy pouch, closed, for use on barrier with locking flange, with filter (two-piece), each	60 ea/mo
A4424	Ostomy pouch, drainable, with barrier attached, with filter (one-piece), each	20 ea/mo
A4425	Ostomy pouch, drainable, for use on barrier with non-locking flange, with filter (two-piece system), each	20 ea/mo
A4426	Ostomy pouch, drainable, for use on barrier with locking flange (two-piece system), each	20 ea/mo
A4427	Ostomy pouch, drainable, for use on barrier with locking flange, with filter (two-piece system) each	20 ea/mo
A4428	Ostomy pouch, urinary, with extended wear barrier attached, with faucet-type tap with valve (one-piece), each	*
A4429	Ostomy pouch, urinary, with barrier attached, with built-in convexity, with faucet-type tap with valve (one-piece), each	20 ea/mo
A4430	Ostomy pouch, urinary, with extended wear barrier attached, with built-in convexity, with faucet-type tap with valve (one-piece), each	*
A4431	Ostomy pouch, urinary, with barrier attached, with faucet-type tap, with valve (one-piece), each	20 ea/mo
A4432	Ostomy pouch, urinary, for use on barrier with non-locking flange, with faucet-type tap with valve (two-piece), each	20 ea/mo
A4433	Ostomy pouch, urinary, for use on barrier with locking flange (two-piece), each	20 ea/mo
A4434	Ostomy pouch, urinary, for use on barrier with locking flange, with faucet-type tap with valve (two-piece), each	20 ea/mo
A4435	Ostomy pouch, drainable, high output, with extended wear barrier (one-piece system), with or without filter, each	*
A4450	Tape, non-waterproof, per 18 square inches	40 ea (18 sq in)/mo
A4452	Tape, waterproof, per 18 square inches	40 ea (18 sq in)/mo
A4455	Adhesive remover or solvent (for tape, cement, or other adhesive), per ounce	16 oz/6 mo
A4456	Adhesive remover, wipes, any type, each	*
A5051	Ostomy pouch, closed, with barrier attached (one-piece), each	60 ea/mo
A5052	Ostomy pouch, closed, without barrier attached, (one-piece), each	60 ea/mo
A5053	Ostomy pouch, closed, for use on faceplate, each	60 ea/mo
A5054	Ostomy pouch, closed, for use on barrier with flange (two-piece), each	60 ea/mo
A5055	Stoma Cap	31 ea/mo
A5056	Ostomy pouch, drainable, with extended wear barrier attached, with filter, (one-piece), each	40 ea/mo
A5057	Ostomy pouch, drainable, with extended wear barrier attached, with built-in convexity, with filter, (one-piece), each	40 ea/mo
A5061	Ostomy pouch, drainable, with barrier attached (one-piece), each	20 ea/mo
A5062	Ostomy pouch, drainable, without barrier attached (one-piece), each	20 ea/mo
A5063	Ostomy pouch, drainable, for use on barrier with flange (two-piece system), each	20 ea/mo
A5071	Ostomy pouch, urinary, with barrier attached (one-piece), each	20 ea/mo
A5072	Ostomy pouch, urinary, without barrier attached (one-piece), each	20 ea/mo
A5073	Ostomy pouch, urinary, for use on barrier with flange (two-piece), each	20 ea/mo

* No established usual maximum, refer to Medicare Utilization Guidelines.

HCPSC	Long Description	Quantity
A5081	Continent device, plug for continent stoma	31 ea/mo
A5082	Continent device, catheter for continent stoma	1 ea/mo
A5083	Continent device, stoma absorptive cover for continent stoma	150 ea/mo
A5093	Ostomy accessory, convex insert	10 ea/mo
A5102	Bedside drainage bottle with or without tubing, rigid or expandable, each	2 ea/6 mo
A5120	Skin barrier, wipes or swabs, each	150 ea/6 mo
A5121	Skin barrier, solid, 6 x 6 or equivalent, each	20 ea/mo
A5122	Skin barrier, solid, 8 x 8 or equivalent, each	20 ea/mo
A5126	Adhesive or non-adhesive, disk or foam pad	20 ea/mo
A5131	Appliance cleaner, incontinence and ostomy appliances, per 16 ounce	1 ea (16 oz btl)/mo
A9270	Non-covered item or service	*

*No established usual maximum, refer to Medicare Utilization Guidelines.



DID YOU KNOW...

Supplier Matching

“Having a supplier that can provide your preferred products when you need them is important. If you are having a difficult time finding a supplier who provides your preferred products, please call 888-808-7456 to learn your options.”


– Secure Start Member

Hollister Secure Start services can connect patients with a supplier of their choice based on their insurance plan for the products you recommend. The amount of supplies they need may vary from the maximum allowable quantities, but this is determined by a physician and should be included on the patient's order and medical record. For additional assistance, visit www.securestartservices.com or call us at 1.888.808.7456.



Index • HCPCS Codes

<u>HCPCS</u>	<u>Page(s)</u>	<u>HCPCS</u>	<u>Page(s)</u>
A4357	83	A5061	50, 51, 71, 73
A4362	73, 80, 83	A5062	73
A4363	82	A5063	21, 38, 39, 43, 44, 45, 75
A4364	81	A5071	74
A4367	76, 80, 87	A5120	81
A4368	83, 88	A5122	81
A4371	80, 83	A5131	82
A4373	33, 35	A6154	63
A4385	23, 78, 80, 81	A4358	84
A4387	89	A5114	85
A4388	17, 48, 49, 52, 71	A4331	85
A4389	88, 89	A9270	84
A4390	19, 56, 58, 59, 61, 62		
A4394	80, 82		
A4397	86, 87		
A4398	86		
A4399	86		
A4402	87		
A4405	83		
A4406	80		
A4407	15, 16, 32, 34, 35, 36		
A4408	15, 32		
A4409	15, 21, 29, 30, 31, 43, 44, 46, 86		
A4410	15, 20		
A4411	23, 78		
A4412	40		
A4414	21, 29, 30, 31, 43, 45, 75		
A4415	29, 30, 43, 45		
A4416	64, 65		
A4417	89		
A4419	40, 41, 75		
A4421	82, 83		
A4424	49, 50, 51, 53, 74		
A4425	38, 39, 75		
A4428	66, 67, 71		
A4429	89		
A4430	20, 68, 69, 70		
A4432	21, 42, 44, 46, 76		
A4435	63		
A4455	81		
A4456	81		
A5052	73		
A5054	41		
A5055	65, 87		
A5056	17, 48, 51		
A5057	18, 19, 55, 56, 57, 69		



The spirit of Hollister
remains the same
today as it did in 1921.

Our Mission

To make life more rewarding
and dignified for the people
who use our products
and services.

Our Vision

To grow and prosper as an
independent, employee-owned
company and in the process
become better human beings.

Our Heritage

The roots of Hollister Incorporated date
back to 1921, when 23-year-old entrepreneur
John Dickinson Schneider founded a small
printing company in Chicago, which he named
JDS Printer Craftsman. From the very beginning,
John Schneider's business philosophy was "only
first class is good enough," and the company
became known for its quality products, quality
service, and quality people.

The Hollister logo, Adapt, Adapt CeraRing, AF300, CeraPlus, ComfortWear, FlexTend, FlexTend M, FlexWear, FormaFlex, HolliHesive, Karaya 5, Lock 'n Roll, m9, New Image, Pouchkins, Premier, Premium, QuietWear, Secure Start, SoftFlex, and "Healthy skin. Positive outcomes." are trademarks of Hollister Incorporated.

Remois is a technology of Alcare Co., Ltd.

All other trademarks and copyrights are the property of their respective owners.

Contact Us

Monday - Friday

8:00am - 5:00pm (CT)

.....
For consumer or healthcare professional questions:

1.888.808.7456

.....
For dealer, distributor, institution and
pharmacy questions or orders:

1.800.323.4060

FAX 1.847.680.1017

.....

Hollister Secure Start Services

Enroll on the Go Clinician Assisted Enrollment

Access to personalized, no-cost support for ostomy patients is easier—with mobile device enrollment! Call us with any questions at 1.800.808.7456.



Hollister Incorporated
2000 Hollister Drive
Libertyville, Illinois 60048
1.800.323.4060



www.hollister.com