




# Advance Plus

Intermittent Catheter System

## Touch Free Protective Tip

- Helps protect the catheter tip from exposure to bacteria located within the first 15mm of the distal urethra
- Helps reduce the risk of carrying bacteria up into the urinary system
- Serves as a guide for proper and easy insertion



## User-Regulated Gel Reservoir

- No additional lubrication required
- Stabilizes catheter for easy insertion and withdrawal

## Additional Trusted Features

- Not made with natural rubber latex
- Environmentally responsible

## Advance Plus Pocket™ Intermittent Catheter System



- Designed for confidence and convenience
- The compact packaging helps make the catheter easier to carry and store
- The same detailed features as the Advance Plus intermittent catheter

## Hollister Incorporated

2000 Hollister Drive, Libertyville, IL USA

[www.hollister.com](http://www.hollister.com)



## Advance Plus Intermittent Catheter System

Advance Plus Pocket	Advance Plus	Advance Plus Kit <sup>†</sup>	Size	Description
–	94064	96064	6 Fr	Firm
–	94084	96084	8 Fr	Firm
93104	94104	96104	10 Fr	Firm
32124	94124	96124	12 Fr	Firm
93144	94144	96144	14 Fr	Firm
93164	94164	96164	16 Fr	Firm
–	94184	96184	18 Fr	Firm
–	95124	97124	12 Fr	Coudé
–	95144	97144	14 Fr	Coudé
–	95164	97164	16 Fr	Coudé

<sup>†</sup>Kit Includes: Advance Plus intermittent catheter system, gloves, underpad drape, and Bensalkonium swabsticks

Qty: 100/box

For detailed clinical questions concerning our products:

**1.888.740.8999**

For orders only: **1.800.323.4060**

For more information or to try a free sample, please visit us at [www.hollister.com/advance](http://www.hollister.com/advance)

