

Maintaining the  
Right Environment  
**Makes all the  
Difference**



The **CeraPlus** Skin Barrier  
*with Remois Technology\**

Available on **Premier** One-Piece and  
**New Image** Two-Piece Pouching Systems



## Peristomal Skin Deserves Better

The importance of peristomal skin health in ostomy care can't be overstated. The condition of this skin greatly affects how well the skin barrier can adhere, as well as the general comfort and well-being of a patient living with a stoma.

Unfortunately, peristomal skin problems are too common. About half of those living with a stoma will report peristomal skin issues at some point in their lives.

Many even think that skin problems are a normal part of living with a stoma.

We disagree. Patients do not need to suffer, and our goal is to help maintain healthy peristomal skin from the start.

## The CeraPlus Skin Barrier

*Available on Premier One-Piece and  
New Image Two-Piece Pouching Systems*

The CeraPlus skin barrier is infused with ceramide, a naturally occurring lipid in the skin, which protects against dryness, and features a formulation to help:


- Protect the skin's natural moisture barrier
- Maintain good peristomal skin health
- Decrease transepidermal water loss (TEWL) from damaged or eroded skin



# The Ceramide Advantage

Ceramide is a natural component of human skin, consisting of waxy lipid molecules. Ceramide links cells of the outer skin (epidermis) together to form a waterproof, protective barrier.

Ceramides help prevent water loss that can lead to skin damage and dryness. Many skin moisturizing products contain ceramide as a component to help protect the skin against loss of its natural moisture.



## Ceramide:

- Is a natural component of the human skin
- Links the cells of the stratum corneum together to form a waterproof, protective barrier
- Helps prevent water loss that can lead to skin damage and dryness

---

*82% of responding clinicians reported they were likely to continue to use and recommend the CeraPlus skin barrier for the patient.\*\**

---

Our new CeraPlus skin barrier is designed to support healthy skin from the start—so your patients maintain good skin health long after they leave your care.

Experience the CeraPlus skin barrier difference.

Contact your Hollister representative for details and a trial sample.

\*\* Results from a 2015 product evaluation (N=193) of the CeraPlus skin barrier in The United States.

# Premier CeraPlus One-Piece Pouching System



## Premier CeraPlus One-Piece Pouching Systems

- Drainable Pouches
- With Tape Border
- Lock 'n Roll Microseal Closure
- Beige with Viewing Option - **With Integrated AF300 Filter**
- Ultra-Clear - **Without Integrated Filter**

### Flat Skin Barriers (Cut-to-Fit) Box of 10

Ultra-Clear Stock No <small>HCPCS: A4388</small>	Beige w/Viewing Option <small>HCPCS: A5056</small>	For Stomas Up To	Pouch Length
8931	8901	2 ½" (64 mm)	12" (30 cm)

### Flat Skin Barriers (Pre-Sized) Box of 10

Beige w/Viewing Option <small>HCPCS: A5056</small>	Barrier Opening	Pouch Length
8925	1" (25 mm)	12" (30 cm)
8930	1 ⅜" (30 mm)	12" (30 cm)
8935	1 ⅝" (35 mm)	12" (30 cm)
8940	1 ⅞" (40 mm)	12" (30 cm)

### Convex Skin Barriers (Cut-to-Fit) Box of 5

Ultra-Clear Stock No <small>HCPCS: A4390</small>	Beige w/Viewing Option <small>HCPCS: A5057</small>	For Stomas Up To	Pouch Length
8984	8914	1" (25 mm)	12" (30 cm)
8988	8918	1 ½" (38 mm)	12" (30 cm)
89811	89111	2" (51 mm)	12" (30 cm)

### Convex Skin Barriers (Pre-Sized) Box of 5

Beige w/Viewing Option <small>HCPCS: A5057</small>	Barrier Opening	Pouch Length
8990	¾" (19 mm)	12" (30 cm)
8991	⅞" (22 mm)	12" (30 cm)
8992	1" (25 mm)	12" (30 cm)
8993	1 ⅛" (29 mm)	12" (30 cm)

### Soft Convex Skin Barriers (Cut-to-Fit) Box of 5

Beige w/Viewing Option <small>HCPCS: A5057</small>	For Stomas Up To	Pouch Length
8954	1" (25 mm)	12" (30 cm)
8958	1 ½" (38 mm)	12" (30 cm)
89511	2 ⅛" (55 mm)	12" (30 cm)

### Soft Convex Skin Barriers (Pre-Sized) Box of 5




Beige w/Viewing Option <small>HCPCS: A5057</small>	Barrier Opening	Pouch Length
8960	¾" (19 mm)	12" (30 cm)
8961	⅞" (22 mm)	12" (30 cm)
8962	1" (25 mm)	12" (30 cm)
8963	1 ⅛" (29 mm)	12" (30 cm)

# New Image CeraPlus Two-Piece Skin Barriers






## New Image CeraPlus Two-Piece Barriers (Without Tape Border) **Box of 5**

### Flat Skin Barriers (Cut-to-Fit)




Stock No	HCPCS 15102/15103: A4409 HCPCS 15104: A4410	Color Match	Flange Size	For Stomas Up To
15102			1¾" (44 mm)	1¼" (32 mm)
15103			2¼" (57 mm)	1¾" (44 mm)
15104			2¾" (70 mm)	2¼" (57 mm)

### Convex Skin Barriers (Cut-to-Fit)




Stock No	HCPCS 15302/15303: A4407 HCPCS 15304: A4408	Color Match	Flange Size	For Stomas Up To
15302			1¾" (44 mm)	1" (25 mm)
15303			2¼" (57 mm)	1½" (38 mm)
15304			2¾" (70 mm)	2" (51 mm)

## New Image CeraPlus Two-Piece Barriers (With Tape Border) **Box of 5**





### Flat Skin Barriers (Cut-to-Fit)

Stock No	HCPCS: A4409	Color Match	Flange Size	For Stomas Up To
11202			1¾" (44 mm)	1¼" (32 mm)
11203			2¼" (57 mm)	1¾" (44 mm)
11204			2¾" (70 mm)	2¼" (57 mm)

### Convex Skin Barriers (Cut-to-Fit)

Stock No	HCPCS: A4407	Color Match	Flange Size	For Stomas Up To
11402			1¾" (44 mm)	1" (25 mm)
11403			2¼" (57 mm)	1½" (38 mm)
11404			2¾" (70 mm)	2" (51 mm)

### Convex Skin Barriers (Pre-Sized)

Stock No	HCPCS: A4407	Color Match	Flange Size	Barrier Opening
11503			1¾" (44 mm)	7⁄8" (22 cm)
11504			1¾" (44 mm)	1" (25 mm)
11505			2¼" (57 mm)	1½" (29 cm)
11506			2¼" (57 mm)	1¼" (32 mm)

Hollister Ostomy. **Details Matter.**

**Hollister Incorporated**  
2000 Hollister Drive  
Libertyville, Illinois 60048  
1.800.323.4060

[www.hollister.com](http://www.hollister.com)



The Hollister logo, AF300, CeraPlus, Lock 'n Roll, New Image, Premier and "Hollister Ostomy. Details Matter." are trademarks of Hollister Incorporated.

\*Remois is a technology of Alcare Co., Ltd.

© 2016 Hollister Incorporated.